Lumax 540

ICD Family

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ICD-familie

ICD Familie

Familia DAI

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System Description

Intended Medical Use

Lumax is the name of a family of implantable cardioverter-defibrillators (ICDs).

Primary objective of the therapy is to prevent sudden cardiac death. The aim is to automatically detect and terminate cardiac arrest caused by ventricular tachyarrhythmia. All major therapeutic approaches from the field of cardiology and electrophysiology are contained within the Lumax family.

Furthermore, the device is capable of treating bradycardia arrhythmias and congestive heart failure with multisite ventricular pacing – known as cardiac resynchronization therapy.

The integrated Home Monitoring component can provide information about occurring rhythm disturbances and delivered therapies close to real time as well as by IEGM Online HD®. Furthermore, statistical data about the patient's condition as well as information about the integrity status of the device itself are sent.

The implantation of an ICD is a symptomatic therapy with the following objectives:

- Termination of spontaneous ventricular fibrillation (VF) through shock delivery
- Termination of spontaneous ventricular tachycardia (VT) through antitachycardia pacing (ATP); in case of ineffective ATP or hemodynamically not tolerated VT, with shock delivery
- Cardiac resynchronization through multisite ventricular pacing (triple-chamber device)
- Compensation of bradycardia through ventricular (single-chamber device) or AV sequential pacing (dual- and triple-chamber device)

Required expertise

In addition to having basic medical knowledge, the user must be thoroughly familiar with the operation of a device system. Only qualified medical specialists having the special knowledge required for the proper use of devices are permitted to use them. If users do not possess this knowledge, they must be trained accordingly.

ICD System

Lumax

The Lumax ICD system consists of the following:

- Single-, dual- or triple-chamber device with connections for bipolar sensing and pacing as well as connections for shock delivery
- ICD leads:
 - One bipolar ICD lead with one or two shock coils for the ventricle (single-chamber device)
 - One bipolar lead for the atrium and one bipolar ICD lead for the ventricle respectively with one or two shock coils (dual-chamber or triple-chamber device)
 - One uni- or bipolar CS lead (coronary sinus lead for the triple-chamber device)
- Programmer with current device program

Note: A tripolar or quadrupolar RV lead has bipolar electrodes as well as 1 or 2 shock coils.

Device

The device's housing is made of biocompatible titanium, welded from outside and thus hermetically sealed. It serves as a potential antipole during shock delivery. The ellipsoid shape facilitates implantation in the pectoral muscle area.

The labeling provides information about the device type and the arrangement of the connections.

Lumax family

The following types with Home Monitoring are available (some device types are not available in all countries):

Device	High energy type: max. 40 J
Single-chamber	Lumax 540 VR-T Lumax 540 VR-T DX
Dual-chamber	Lumax 540 DR-T
Triple-chamber	Lumax 540 HF-T

NBD and NBG codes

VVE is the NBD code for the antitachycardia mode of the single-chamber, dual-chamber and triple-chamber devices:

V	Shock in the ventricle
V	Antitachycardia pacing (ATP) in the ventricle
E	Detection via IEGM analysis

DDDR is the NBG code for the antibradycardia mode of the dual-chamber and triple-chamber devices:

D	Pacing in both chambers
D	Sensing in both chambers
D	Pulse inhibition and pulse triggering
R	Rate adaptation

VWIR is the NBG code for the antibradycardia modes of the single-chamber devices:

V	Ventricular pacing
V	Sensing in the ventricle
I	Pulse inhibition in the ventricle
R	Rate adaptation

ICD leads

The leads are sheathed with biocompatible silicone. They can be flexibly maneuvered, are long-term stable, and are equipped for active or passive fixation. They are implanted using a lead delivery set. Some leads are coated with polyurethane which is known to increase the sliding properties for the lead. Leads with steroids reduce inflammatory processes. The fractal design of the electrodes provides for low pacing thresholds.


Note: Device and lead have to match; Kentrox A+ Steroid and Linox Smart DX leads match the Lumax 540 VR-T DX.

Device/ICD lead connection scheme

BIOTRONIK ICDs are designed for connection to ICD leads with a bipolar IS-1 connection (pacing/sensing) and DF-1 connection for shock coils.

Note: A bipolar as well as a unipolar coronary sinus lead can be connected to the left ventricular CS connection of the Lumax HF variant. Several pacing polarities can be set with a bipolar CS lead.

Connection schemes of the device types:

Device and lead connections:

	VR-T	VR-T DX	DR-T	HF-T
Atrium	—	IS-1 bipolar		
Right ventricle	IS-1 bipolar 2 x DF-1 unipolar			
Left ventricle	—		IS-1 unipolar or bipolar	

Possible technical failures

In principle, technical failures in the device due to component failures cannot be ruled out; they occur, however, very rarely. Possible causes of such malfunctions include the following:

- Battery depletion
- Lead dislodgment
- Lead fracture
- Insulation defects

Programmer

The portable programmer is used to transfer the current device program to the device. In addition to this, the programmer is used for interrogation and storage of data from the device. The programmer acts as an ECG and IEGM monitor with Miniclinic.

The programmer communicates with the device via the programming head. It has a TFT touch screen with a color display, on which the ECG, IEGM, marker and functions are shown simultaneously.

The programmer has, among others, the following functions:

- Perform all tests during in-office follow-up
- Display and print real-time and saved IEGMs with annotated markers
- Determine the pacing threshold

BIOTRONIK Home Monitoring®

In addition to effective pacing therapy, BIOTRONIK provides a complete therapy management system:

- With Home Monitoring, diagnostic and therapeutic information and technical data are automatically sent to a stationary or mobile transmitter via an antenna in the device header. The data are encrypted and sent from the transmitter to the BIOTRONIK Service Center via the cellular phone network.
- The received data are deciphered and evaluated. Each physician can set the criteria for evaluation to be used for each patient and can configure the time of notification via E-mail, SMS or fax.
- A clear overview of the analysis results is displayed for the attending physicians on the protected Internet platform HMSC (Home Monitoring Service Center).
- Data transmission from the device is performed with a daily device message.
- Device messages, which indicate special events in the patient's heart or in the device, are forwarded immediately.

Technical manuals

The following technical manuals provide information about usage of the device system:

- Technical manual for the device
- Technical manual for the HMSC
- Technical manual for the programmer
- User manual for the device program:
 - As a help function in the user interface
 - As a file on CD
- Technical manual for the leads
- Technical manual for cables, adapters and accessories

Therapeutic and Diagnostic Functions

Overview

Depending on the type, the device program contains not only the ICD functions but also all PM functions for 1, 2, or 3 chambers. The patient's heart rhythm is continuously monitored; each arrhythmia is classified according to the heart rate and the adjustable detection criteria. Depending on the preset values, antibradycardia as well as anti-tachycardia therapy is inhibited or delivered.

Antitachycardia pacing

The ICD can treat ventricular tachycardia with antitachycardia pacing (ATP); ATP can also be delivered in the VF zone (ATP One Shot) when the rate stability indicates that this will be effective before shock delivery (monomorphic rapid VTs).

Cardioversion, defibrillation

The ICD can treat ventricular tachyarrhythmia with cardioversion and/or defibrillation. Shock polarity and energy can be programmed for individual therapies. Shock energies between 1.0 and 40 J are possible. Before delivery of the first shock, the ICD can verify that the tachyarrhythmia is still ongoing; in this verification mode the device can identify spontaneous conversion of the tachyarrhythmia and cancel shock delivery.

- The shock paths can be set between the different shock coils (SVC/RV) and/or the housing.

Home Monitoring: Message retrieval

The device sends information to the service center in two ways; automatically once a day, and when triggered by certain events. All data can be viewed on the Internet subsequently.

- IEGM-Online HD® with 3 channels in high resolution (High Density) with annotated markers
- Direct information on sustained atrial arrhythmia including IEGM Online HD
- Regular IEGM recordings of the intrinsic rhythm for remote diagnosis
- Daily transmission of atrial heart rate variability (SDANN)

Transmission of the following data:

- Device settings
- Impedance measurement values
- Sensing measurement values
- 30-second IEGM recording
- Statistics

Transmission is carried out on a weekly or monthly basis.

Antibradycardia pacing

Innovative rate hystereses, automatic sensor functions, and a night program promote the patient's intrinsic rhythm, avoid overdrive pacing, and facilitate adaptation of the device to the individual needs of the patient.

Positive AV hysteresis functions support the intrinsic conduction and thus the natural contraction sequence. Negative AV hysteresis functions support the cardiac resynchronization therapy by maintaining pacing in stress situations.

- For resynchronization of the ventricles, 3-chamber devices have functions for multisite ventricular pacing with possible WV delays in either direction.
- To ensure that no additional surgery is necessary in case of a left-sided increase of pacing threshold or undesired phrenic nerve stimulation by the CS lead, different pacing polarities can be set for the CS lead with a triple-chamber device.
- Setting an upper pacing rate for the atrium prevents unspecific atrial pacing, thus improving the termination of pacemaker-mediated tachycardia.
- Automatic threshold monitoring (ATM) for the left and right ventricle

Diagnostic functions

The device can record several types of arrhythmia episodes and store them in the Holter memory:

- VT/VF therapy episodes
- SVT episodes (requirement SMART)
- VT monitoring zones
- AT/AF monitoring zones

The device stores diagnostic data about arrhythmia episodes:

- Detection counter and therapy counter
- Temporal characteristics of the episode
- Details on the last ATP and shock deliveries
- Shock data
- Therapy history
- Triple-channel IEGMs including marker channels, up to 30 min

Comprehensive memory functions such as histograms, Holter, mode switching long- and short-term trends and the activity reports facilitate the evaluation of both patient and device status.

The follow-up assistant functions have been largely automated.

The shock impedance is measured in an automatic procedure involving below-threshold electrical pulses. This process is painless for the patient. The impedance values are automatically measured 4 times a day. They are averaged and provided in a trend in the device and via Home Monitoring.

Painless measurement of shock impedance can also be carried out manually with the programmer.

- Lumax 540 VR-T DX is new (DX stands for diagnostics): Only one lead is required for improved detection of AT and AF even with a single-chamber device. This enables sensing in the atrium and ventricle as well as AV-sequential pacing. The SMART algorithm helps distinguish more clearly between SVT and VT. The IEGM is recorded in both chambers.

Prior to Implantation

Indications

Lumax ICDs can treat life-threatening ventricular arrhythmias with antitachycardia stimulation and defibrillation.

Single- and dual-chamber ICD

Lumax single and dual chamber ICDs are indicated for patients with the following risk:

- Sudden cardiac death caused by ventricular arrhythmias

Triple-chamber ICD

Lumax triple-chamber ICDs are indicated for patients with the following threats:

- Sudden cardiac death caused by ventricular arrhythmias
- Congestive heart failure with ventricular asynchrony

Lumax ICDs are also indicated for primary prophylaxis in congestive heart failure patients.

Contraindications

Known contraindications:

- Tachyarrhythmia caused by temporary or reversible irritation e.g. poisoning, electrolyte imbalance, hypoxia, sepsis or acute myocardial infarction
- Such frequent VT or VF that the therapies would cause an unacceptably rapid depletion of the device batteries
- VT with few or without clinically relevant symptoms
- VT or VF treatable by surgery
- Concomitant diseases that would substantially limit a positive prognosis
- Accelerated idioventricular rhythm

Shipping and Storage

Information on the package

The storage package for the device comprises a folded cardboard box with a quality control seal and a label.

- Device name with serial number and PID number
- Technical details
- Details on sterility and use by date
- Ambient conditions during shipping and storage

Ambient conditions

5°C to 45°C

Storage location

- Devices are not to be stored close to magnets or sources of electromagnetic interference.

Implantation

Sterility

Delivery

The device and the accessories have been gas sterilized. Sterility is guaranteed only if the plastic container and quality control seal have not been damaged.

Sterile container

The device and accessories are packaged respectively in two separately sealed plastic containers. The inner plastic container is also sterile on the outside so that it can be transferred in a sterile state during implantation.

Preparing the Implantation

Having the required parts ready

Note: Ensure that sterile spare parts are available for all parts that are to be implanted.

- Device with blind plug DF-1 and IS-1
- BIOTRONIK ICD leads complying with the requirements of EC Directive 90/385/EEC, and lead delivery set
- BIOTRONIK programmer complying with the requirements of EC Directive 90/385/EEC
 - ECG cable PK-222
 - 2 PK-141 patient cables or 2 PK-67 patient cables, either with patient adapters PA-2 or PA-4 for connection of the pacing and sensing leads (IS-1) to ICS (dual- and triple-chamber devices)
 - 1 PK-144 patient cable with patient adapter PA-3 for connection of shock coils (DF-1) to the programmer
- External multi-channel ECG recorder/monitor
- External defibrillator and paddles or adhesive electrodes.

WARNING

Non-terminated ventricular arrhythmia

If the ICD does not deliver an adequate therapy, ventricular arrhythmia is not terminated.

- Keep an external defibrillator ready.

Unpacking the device

Proceed as follows:

WARNING

Inadequate therapy due to previously damaged device

If the unpacked device is dropped on a hard surface during handling, discontinue use and return it to BIOTRONIK. Use a replacement device.

- 1 Peel off the sealing paper of the non-sterile outer plastic container at the marked position in the direction indicated by the arrow.
The inner plastic container may not come into contact with persons who have not sterilized their hands or gloves, nor with non-sterile instruments.
- 2 Take hold of the inner plastic container by the gripping tab and take it out the outer plastic container.
- 3 Peel off the sealing paper of the sterile inner plastic container at the marked position in the direction indicated by the arrow.

Connecting ICD Leads

Precautionary measures

Note: Use only adapters approved by BIOTRONIK for ICD leads with different connections. Should you have questions about the compatibility of other manufacturers' ICD leads, please contact BIOTRONIK.

WARNING

Short circuit due to open lead connections

Open – and thus not electrolyte-proof – IS-1 or DF-1 connections may cause undesired current flows to the body and penetration of body fluid into the device.

- Close unused DF-1 connections with the blind plug DF-1, unused IS-1 connections with the blind plug IS-1.

⚠ WARNING

Far-field sensing or insufficient defibrillation

When shock coils and pace/sense electrodes are not spaced sufficiently apart or positioned inappropriately, this can lead to far-field sensing or insufficient defibrillation:

- The distance between two shock coils must be greater than 6 cm.
- Pace and sense electrodes must not contact each other.

⚠ CAUTION

Damage to the header while handling the blind plug

For every connection there is a blind plug in the header; the provided set screws must be carefully loosened or tightened.

- Loosen set screws with the supplied screwdriver. Use BIOTRONIK screwdrivers with torque control only!
- Do not forcibly pull out the blind plug!
- If lead repositioning is necessary, re-order sterile screwdrivers from BIOTRONIK.

⚠ WARNING

Malfunctioning of the therapy functions or triggering of tachycardia by older, unused leads

If, upon replacing the device, predecessor leads are no longer used but left where they are, then their connections have to be insulated so that there is no additional uncontrolled current path to the heart.

Connecting the ICD lead connector to the device

Leads are connected to the device according to the diagram on the header. Proceed as follows:

- 1 Disconnect styles and insertion aids from the lead connector.
- 2
 - Connect ventricle shock coil to RV.
 - Connect vena-cava shock coil or subcutaneous shock coil to SVC.
- 3
 - Connect the bipolar IS-1 lead connector atrium to A S or P/S A.
 - Connect the bipolar IS-1 lead connector ventricle to P/S V.
 - Connect the unipolar or the bipolar IS-1 lead connector of the CS lead to P/S CS.
- 4 Push the lead connector into the header without bending the conductor until the connector tip becomes visible behind the set screw block.

- 5 If the lead connector cannot be inserted completely, the set screw may be protruding into the cavity of the set screw block.
Carefully loosen the set screw without completely unscrewing it, so that it does not become tilted upon retightening.
- 6 Use the screwdriver to perpendicularly pierce through the slitting in the center of the silicone plug until it reaches the set screw.
- 7 Turn the set screw clockwise until the torque control starts [you will hear a clicking sound].
- 8 Carefully withdraw the screwdriver without retracting the set screw.
 - Bipolar IS-1 connectors: tighten both set screws!
 - When you withdraw the screwdriver, the silicone plug automatically seals the lead connection safely.

Setting the lead polarity manually

Due to the risk of an entrance/exit block, bipolar lead polarity [sensing/pacing] should only be set if bipolar leads are implanted.

Implanting the ICD

Implantation site

In general the ICD is implanted subpectorally on the left depending on the lead configuration as well as the anatomy of the patient.

ICD status prior to implantation

The ICD is shipped in a deactivated mode and can be connected and implanted in this state.

⚠ CAUTION

Shipment mode active!

During shipment, the device has been set to shipment mode. The shipment mode reduces the battery's energy consumption, which may increase longevity. The shipment mode controls the charge time of automatic capacitor reformations.

The shipment mode must be deactivated before the conclusion of the implantation procedure.

The shipment mode will be automatically deactivated as soon as electrophysiological tests (e.g. impedance measurement) have been performed.

⚠ WARNING

Unintentional shock delivery

A shock could be delivered when handling an activated ICD.

- If an already activated ICD is to be implanted: turn off ICD therapy.

⚠ WARNING

Interference of heart rate monitoring

Radio frequency as during electrocautery can interfere with electromyographic monitoring, e.g. be interpreted as arrhythmia and trigger a shock.

- Turn off the detection function of the ICD during cauterization; the pacemaker function can remain active.
- Additionally check the peripheral pulse of the patient.

Notes on applying a programming head (PGH)

The programming head (PGH) features a diagram of the device to assist in positioning the head.

- Make sure the PGH is positioned correctly to ensure proper telemetry.

Notes on applying a permanent magnet

Applying a permanent magnet interrupts sensing and therapy of tachycardia events.

After 8 hours of this type of deactivation, the device automatically activates the therapy functions again to prevent unintentional permanent deactivation.

- If detection interruptions of longer than 8 hours are required, the magnet has to be briefly removed from the device. The 8 hour countdown restarts when the magnet is applied again.
- Use BIOTRONIK accessories, permanent magnets of type M-50 or programming head PGH with magnet.

Sequence

Proceed as follows:

1	Shape the device pocket and prepare the vein.
2	Implant the ICD leads and perform measurements.
3	Connect the device and ICD leads.
4	Insert the device.
5	Check the device with standard tests.
6	Close the device pocket.

⚠ WARNING

Telemetry interference with inadequate pulse amplitude

If the pulse amplitude of the temporary program is not effective, the patient may get into a hemodynamically critical condition. Telemetry interference can prevent correction of the critical amplitude.

- Lift the programming head 30 cm; the device switches automatically to the permanent program.

⚠ WARNING

Loss of pacing during exclusively left ventricular pacing

When only LV pacing is set and lead displacement occurs, the following dangers exist: loss of ventricular pacing and the ATP therapy as well as induction of atrial arrhythmia.

- Carefully check the pacing parameters.

Activating the ICD

Proceed as follows:

1	Select the device program for Lumax on the programmer.
⚠ WARNING	
2	Activate ICD therapy.

After Implantation

Follow-up

Follow-up intervals

Follow-ups must be performed at regular, agreed intervals.

- The first follow-up should be carried out by the physician using the programmer (in-office follow-up) approximately 3 months after implantation following the lead ingrowth phase.
- The next in-office follow-up should be carried out once a year and no later than 12 months after the last in-office follow-up.

Follow-up with BIOTRONIK Home Monitoring®

Monitoring using the Home Monitoring function does not serve to replace regular in-office appointments with the physician required for other medical reasons. Follow-up supported by Home Monitoring can be used to functionally replace in-office follow-up under the following conditions:

- The patient was informed that the physician must be contacted if symptoms worsen or if new symptoms arise despite use of the Home Monitoring function.
- Device messages are transmitted regularly.
- The physician decides whether the data transmitted via Home Monitoring with regard to the patient's clinical condition as well as the technical state of the device system are sufficient. If not, an in-office follow-up has to be carried out.

Possible early detection due to information gained via Home Monitoring may necessitate an additional in-office follow-up. For example, the data may indicate at an early stage lead problems or a foreseeable end of service time (ERI). Furthermore the data could provide indications of previously unrecognized arrhythmias or modification of the therapy by reprogramming the device.

Follow-up with the programmer

Use the following procedure for in-office follow-up:

1	Record and evaluate the external ECG.
2	Check the sensing and pacing functions.
3	Interrogate the device.
4	Evaluate the status and automatically measured follow-up data.
5	Possibly evaluate statistics and Holter/IEGM recording.
6	Manually perform standard tests if necessary.
7	Possibly customize program functions and parameters.
8	Transmit the program permanently to the device.
9	Print and document follow-up data (print report).
10	Finish the follow-up for this patient.

Note for physicians

Use the patient ID card to obtain information on magnet behavior and tachyarrhythmia detection.

⚠ WARNING

Possible loss of left ventricular pacing if the pacing threshold has been determined exclusively by ATM

The pacing threshold determined by ATM should not be used directly to determine the left ventricular pacing amplitude (LV). The effectiveness of the LV pacing amplitude has to be confirmed.

Notes for patients

A patient brochure and a patient ID card are supplied with the device.

- Provide the patient with the patient brochure and patient ID card.
- Draw the patient's attention to prohibitory signs:
 places with prohibitory signs must be avoided.

Replacement Indication

Foreword

The battery level can be continuously monitored with Home Monitoring and checked during follow-ups.

Possible battery levels

- BOL: Beginning of Life (i.e. BOS: Beginning of Service): > 70% charge
- MOL 1: Middle of Life : 70% to 40% residual charge
- MOL 2: Middle of Life : < 40% residual charge
- ERI: Elective Replacement Indication, (i.e. RRT: Recommended Replacement Time)
- EOS: End of Service

ERI

The device can monitor the cardiac rhythm for at least 3 months and deliver at least six maximum energy shocks until EOS.

The selected parameters in the device program do not change.

⚠ CAUTION

Temporally limited therapy

When the ERI display is noticed for the first time during a follow-up, the remaining service time can be much less than 3 months even if follow-ups have been performed at three month intervals.

Action on ERI

Proceed as follows:

1	Replace device soon.
---	----------------------

EOS

- End of Service can be detected by Home Monitoring.
- VT and VF detection and all therapies are deactivated.
- The antibradycardia function remains active in the VVI mode:
 - Basic rate 50 ppm
 - Without special pacemaker functions, such as hysteresis, etc.
 - Amplitude of 6 V
 - Pulse width of 1.5 ms

⚠ WARNING

EOS before explantation

If EOS occurs before replacement of the device, the patient is without antiarrhythmic therapy and thus at risk of a life threatening condition!

Action on EOS

Proceed as follows:

1	Replace device immediately.
2	Monitor patient constantly until immediate replacement of the device!

Explantation and Disposal

Note: Before the cremation of a deceased patient, the device must be explanted.

⚠ WARNING

Unintentional shock delivery

A shock could be delivered when handling an activated ICD. Devices may not be explanted when detection is activated.

- Turn off the detection function before explantation.
- Do not simply cut the ICD leads.

Explanting the ICD

Proceed as follows:

1	Interrogate the device status.
2	Deactivate VT and VF therapies.
3	Disconnect the leads from the header.
4	Remove the device and leads using state-of-the-art technology.
5	Dispose of the device in an environmentally safe manner.
6	Fill out explantation form and send to BIOTRONIK.

Note: Normal oxidation processes may cause ICD housing discolorations. This is neither a device defect nor does it influence device functionality.

Disposing of the ICD

Proceed as follows:

1	Clean the explant with an at least 1% sodium-hyperchlorine solution.
2	Rinse off with water.
3	Send to BIOTRONIK for disposal in an environmentally safe manner.

Cautionary Notes

Medical Complications

General information on possible complications

It is impossible to guarantee the efficacy of antitachycardia or fibrillation therapy, even if these therapies have proved successful during intraoperative tests or subsequent electrophysiological tests. In rare circumstances, the set therapies may prove ineffective or, in unfavorable conditions, even life-threatening. It is possible for therapies to induce or accelerate tachycardia and cause sustained ventricular flutter or fibrillation.

Known possible complications:

The following are some of the known medical complications related to pacemaker therapy and ICD therapy:

- Formation of necrotic tissue
- Vascular damage
- Thrombosis

- Embolism
- Increase in pacing threshold
- Detection failures
- Foreign body rejection
- Cardiac tamponade
- Muscle or nerve stimulation
- Device-induced arrhythmias
- Perforation of device pockets
- Infections
- Psychological intolerance or psychological dependency

Actions to prevent device-induced complications

△ CAUTION

Transmission of atrial tachycardia to the ventricle

Set the following parameters to avoid unphysiological pacing in the ventricle in cases of high atrial rates or sinus tachycardia:

- Activate Mode Switching for indicated patients.
- Set the upper rate and the refractory periods to prevent abrupt ventricular rate switching.
- Prefer Wenckebach response and avoid 2:1 behavior.

Set all parameters to prevent constant switching between atrial and ventricle controlled modes to prevent unphysiological rhythm switching upon loss of AV sequential pacing.

△ CAUTION

Pacemaker-mediated tachycardia in cases of retrograde conduction

Pacemaker-mediated tachycardia can occur in patients with retrograde conduction.

- Measure the retrograde conduction time.
- Switch on PMT protection to prevent pacemaker-mediated tachycardia.
- Set the VA criterion.

Skeletal myopotentials

Bipolar sensing and control of sensitivity are adapted by the device to the rate spectrum of intrinsic events so that skeletal myopotentials are usually not detected. Skeletal myopotentials can very rarely be sensed as intrinsic events and, depending on the interference pattern, may provoke inhibition or antiarrhythmia therapy.

Risky Therapeutic and Diagnostic Procedures

External defibrillation

The device is protected against the energy that is normally induced by external defibrillation. Nevertheless, any implanted device may be damaged by external defibrillation. Specifically, the current induced in the implanted leads may result in necrotic tissue formation close to the electrode/tissue interface. As a result, sensing properties and pacing thresholds may change.

Note: Place adhesive electrodes anterior-posterior or perpendicular to the axis formed by the device to the heart at least 10 cm away from the device and from implanted leads.

Due to possible damage to the patient or the device, the use of the following procedures is contraindicated:

- Therapeutic ultrasound and diathermy
 - Damage to the patient via excess warming of body tissue near the device system
- Radiation therapy
 - Shield the device sufficiently against radiation.
 - Check the system for integrity after applying radiation.
 - Radiation can cause latent damage.
- Transcutaneous Electrical Nerve Stimulation (TENS)
- Lithotripsy
- Magnetic resonance imaging and the associated magnetic flux density
 - Damage or destruction of the device system due to strong magnetic interaction
 - Damage to the patient via excess warming of body tissue near the device system
- Electrocautery and high-frequency surgery
 - Damage to the patient via the induction of arrhythmia or ventricular fibrillation
- Hyperbaric oxygen therapy
- Applied pressures higher than normal pressure

△ CAUTION

Interference with the device and risk to the patient due to electrical current during medical treatment

If electrical current from an external source is conducted through the body for diagnostic or therapeutic purposes, then the device has to be switched off or carefully monitored during the initial treatment phases.

Possible Technical Complications

- Device component failures
- Battery depletion
- Lead dislocation
- Lead fracture
- Insulation defect

⚠ CAUTION

Technical failures due to damaged leads

Switch the automatic impedance measurement function on to detect lead errors. Impedance values that indicate a loss of lead integrity are documented in the event list.

Possible Interference

Possible impairments due to EMI

All devices (PM/ICD) are susceptible to interference. The signals generated by sources of interference could be interpreted as heart activity by the device, and/or could disturb measurements required for adapting the rate of the device:

- Depending on the pacing mode and the type of interference, these sources of interference may lead to pacemaker pulse inhibition or triggering, an increase in the sensor-dependent pacing rate, or a fixed-rate pulse delivery.
- Under unfavorable conditions, for example during diagnostic or therapeutic procedures, the interference sources may induce such a high level of energy into the artificial pacing system that the device and/or cardiac tissue around the lead tip is damaged.
- BIOTRONIK devices have been designed so that their susceptibility to EMI is minimized.
- However, due to the variety and intensity of EMI, absolute safety is not possible. It is generally assumed that EMI produces only minor symptoms, if any, in patients.

Note: Indicate the possibility of interference to your patient if this is expected to have a clinically relevant impact. Protect the device against the interference itself or its effect by means of appropriate programming.

Static magnetic fields

Above a flux density of 1.8 mT, the Reed contact in the ICD closes; the ICD therapy is then deactivated. If the magnetic field decreases to below 1 mT, the Reed contact opens and the ICD therapy is active again.

Possible sources of interference

Interference can be caused by:

- Household appliances
- Safety locks or anti-theft installations
- Strong electromagnetic fields
- Cellular phones and transmitters (CardioMessenger)
 - Patients are advised to hold cellular phones to the ear opposite the side on which the device is implanted. Furthermore the cellular phone or transmitter should be held at least 15 cm away.
 - Cellular phones and transmitters emit signals when in stand-by mode, i.e. when not in use. Therefore these devices should not be carried in a breast pocket or within a radius of 15 cm around the implanted device.
 - Electromagnetic interference only has a temporary effect. Implanted devices function again properly when the cellular phone or transmitter is moved away.

Technical Data

Bradycardia Therapy Parameters

Pacing Modes

Lumax 540 family

The following pacing modes are available:

Device type	Pacing modes	Standard
DR-T; HF-T	DDD, DDDR, DDI, DDIR, VDD, VDDR, VDI, VDIR, VI, VIR, AAI, AAIR, OFF	DDD

Device type	Pacing modes	Standard
VR-T	VVI, VVIR, OFF	VVI
VR-T DX	VDD, VDDR, VDI, VDIR, VVI, VVIR OFF	VVI

Timing Parameters Lumax DR-T

Basic rate day/night

Parameter	Range of values	Standard
Basic rate	30 ... [5] ... 100 ... [10] ... 160 ppm	60 ppm
Night rate	OFF 30 ... [5] ... 100 ppm	OFF
Night begins	00:00 ... [1 min] ... 23:59 hh:mm	22:00:00 hh:mm
Night ends	00:00 ... [1 min] ... 23:59 hh:mm	06:00 hh:mm

Rate hystereses

Parameter	Range of values	Standard
Rate hysteresis	OFF -5 ... [-5] ... -90 ppm	OFF
Repetitive hysteresis	OFF 1 ... [1] ... 15	OFF
Scan hysteresis	OFF 1 ... [1] ... 15	OFF

AV delay

Parameter	Range of values	Standard
AV delay	Low; medium; high; fixed; individual	Low
AV delay 1 at rate 1	15 [fixed]; 40 ... [5] ... 350 ms 30 ... [10] ... 150 ppm	- 60 ppm

Parameter	Range of values	Standard
AV delay 2 at rate 2	40 ... [5] ... 350 ms 40 ... [10] ... 160 ppm	- 130 ppm
Sense compensation	OFF -5 ... [5] ... -60 ms	-30 ms
Safety window	100 ms	100 ms

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Parameter	Range of values	Standard
AV hysteresis mode	OFF Positive; Negative; IRSplus (only for DR devices)	OFF
AV hysteresis	10 ... [10] ... 150 ms	50 ms
Positive repetitive AV hysteresis	OFF 1 ... [1] ... 10	OFF
Negative repetitive AV hysteresis	OFF 1 ... [1] ... 15 ... [5] ... 100 ... [10] ... 180	OFF
AV scan hysteresis	OFF 1 ... [1] ... 10 ms	OFF

IRSplus

Parameter	Range of values	Standard
IRSplus	OFF; ON	OFF
AV hysteresis at IRSplus	400 ms	400 ms
Repetitive AV hysteresis at IRSplus	OFF 1 ... [1] ... 10	5
AV scan hysteresis at IRSplus	OFF 1 ... [1] ... 10	5

Post-ventricular atrial refractory period

Parameter	Range of values	Standard
PVARP	AUTO 175 ... [25] ... 600 ms	250 ms

Parameter	Range of values	Standard
PVARP after PVC	PVARP + 225 ms [max: 600 ms] is automatically programmed	475 ms

PVC classification; PVC lock-in protection

Parameter	Range of values	Standard
PVC discrimination after A _s	250 ... [50] ... 450 ms	350 ms

Upper rate (UTR)

Parameter	Range of values	Standard
Upper rate	90 ... [10] ... 160 ppm	130 ppm
Atrial upper rate	OFF 240 ppm	240 ppm

Mode switching

Parameter	Range of values	Standard
Intervention rate	OFF 100 ... [10] ... 250 ppm	160 ppm
Onset criterion X	3 ... [1] ... 8	5
Resolution criterion Y	3 ... [1] ... 8	5
Mode	DDI; DDI(R) when permanent DDD(R) VVI, VVI(R) when permanent VDD(R)	DDI VVI

Post Mode Switching Response (PMSR)

Parameter	Range of values	Standard
Post mode switching rate	OFF +5 ... [5] ... +50 ppm	+10 ppm
Post mode switching duration	1 ... [1] ... 30 min	1 min

PMT protection

Parameter	Range of values	Standard
PMT detection/termination	OFF; ON	ON
VA criterion	250 ... [10] ... 500 ms	350 ms

Pulse amplitude and pulse width

Parameter	Range of values	Standard
Pulse amplitude A	0.2 ... [0.1] ... 6.2; 7.5 V	2.8 V
Pulse width A	0.4; 0.5; 0.7; 1.0; 1.2; 1.5 ms	0.4 ms
Pulse amplitude RV	0.2 ... [0.1] ... 6.2; 7.5 V	2.8 V
Pulse width RV	0.4; 0.5; 0.7; 1.0; 1.2; 1.5 ms	0.4 ms

Automatic threshold monitoring

Automatic threshold monitoring (ATM)

Parameter	Range of values	Standard
ATM RV	ON; OFF	ON

Timing Parameters Lumax 540 VR-T DX

Basic rate day/night

Parameter	Range of values	Standard
Basic rate	30 ... [5] ... 100 ... [10] ... 160 ppm	60 ppm
Night rate	OFF 30 ... [5] ... 100 ppm	OFF
Begin of night	00:00 ... [1 min] ... 23:59 hh:mm	22:00 hh:mm
Night ends	00:00 ... [1 min] ... 23:59 hh:mm	06:00 hh:mm

Rate hystereses

Parameter	Range of values	Standard
Rate hysteresis	OFF -5 ... [-5] ... -90 ppm	OFF
Repetitive hysteresis	OFF 1 ... [1] ... 15	OFF
Scan hysteresis	OFF 1 ... [1] ... 15	OFF

AV delay

Parameter	Range of values	Standard
AV delay	Low; Medium; High; Fixed; Individual	Low
AV delay 1 at rate 1	15 [fixed]; 40 ... [5] ... 350 ms 30 ... [10] ... 150 ppm	- 60 ppm
AV delay 2 at rate 2	40 ... [5] ... 350 ms 40 ... [10] ... 160 ppm	- 130 ppm

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Parameter	Range of values	Standard
AV hysteresis mode	OFF Positive; Negative; IRSplus	OFF
AV hysteresis	10 ... [10] ... 150 ms	50 ms
Positive repetitive AV hysteresis	OFF 1 ... [1] ... 10	OFF
Negative repetitive AV hysteresis	OFF 1 ... [1] ... 15 ... [5] ... 100 ... [10] ... 180	OFF
AV scan hysteresis	OFF 1 ... [1] ... 10	OFF

IRSplus

Parameter	Range of values	Standard
IRSplus	OFF; ON	OFF
AV hysteresis at IRSplus	400 ms	400 ms
Repetitive AV hysteresis at IRSplus	OFF 1 ... [1] ... 10	5
AV scan hysteresis at IRSplus	OFF 1 ... [1] ... 10	5

Post-ventricular atrial refractory period

Parameter	Range of values	Standard
PVARP	AUTO; 175 ... [25] ... 600 ms	250 ms
PVARP after PVC	PVARP + 225 ms (max: 600 ms) is automatically programmed	475 ms

PVC classification; PVC lock-in protection

Parameter	Range of values	Standard
PVC discrimination after A _s	250 ... [50] ... 450 ms	350 ms

Upper rate (UTR)

Parameter	Range of values	Standard
Upper rate	90 ... [10] ... 160 ppm	130 ppm
Atrial upper rate	OFF 240 ppm	240 ppm

Mode switching

Parameter	Range of values	Standard
Intervention rate	OFF 100 ... [10] ... 250 ppm	160 ppm
Activation criterion X	3 ... [1] ... 8	5
Deactivation criterion Y	3 ... [1] ... 8	5
Mode	VDD, VDD(R) when permanent VDD(R)	VDD

Post Mode Switching Response (PMSR)

Parameter	Range of values	Standard
Post mode switching rate	OFF +5 ... [5] ... +50 ppm	+10 ppm
Post mode switching duration	1 ... [1] ... 30 min	1 min

PMT protection

Parameter	Range of values	Standard
PMT detection/termination	OFF; ON	ON
VA criterion	250 ... [10] ... 500 ms	350 ms

Pulse amplitude and pulse width

Parameter	Range of values	Standard
Pulse amplitude RV	0.2 ... [0.1] ... 6.2; 7.5 V	2.8 V
Pulse width RV	0.4; 0.5; 0.7; 1.0; 1.2; 1.5 ms	0.4 ms

Automatic threshold monitoring

Automatic threshold monitoring (ATM)

Parameter	Range of values	Standard
ATM RV	ON; OFF	ON

Timing Parameters Lumax VR-T

Basic rate day/night

Parameter	Range of values	Standard
Basic rate	30 ... [5] ... 100 ... [10] ... 160 ppm	60 ppm
Night rate	OFF 30 ... [5] ... 100 ppm	OFF
Begin of night	00:00 ... [1 min] ... 23:59 hh:mm	22:00 hh:mm
End of night	00:00 ... [1 min] ... 23:59 hh:mm	06:00 hh:mm

Rate hysteresis

Parameter	Range of values	Standard
Rate hysteresis	OFF -5 ... [-5] ... -90 ppm	OFF
Repetitive hysteresis	OFF 1 ... [1] ... 15	OFF
Scan hysteresis	OFF 1 ... [1] ... 15	OFF

Pulse amplitude and pulse width

Parameter	Range of values	Standard
Pulse amplitude V	0.2 ... [0.1] ... 6.2; 7.5 V	2.8 V
Pulse width V	0.4; 0.5; 0.7; 1.0; 1.2; 1.5 ms	0.5 ms

Automatic threshold monitoring

Automatic threshold monitoring (ATM)

Parameter	Range of values	Standard
ATM RV	ON; OFF	OFF

Timing Parameters Lumax HF-T

Ventricular pacing

Parameter	Range of values	Standard
Ventricular pacing	RV; BiV for: DDD(R); DDI(R); VDD(R); VDI(R);VI(R) LV for: DDD(R); VDD(R)	BiV
LVT-wave protection	OFF; ON	ON
Triggering	OFF RVs; RVs+RVES	RVs
Maximum trigger rate	AUTO; 90 ... [10] ... 160 ppm	AUTO

VV delay after Vp

Parameter	Range of values	Standard
Initially paced chamber	RV; LV	LV
VV delay	0 ... [5] ... 100 ms	5 ms
VV delay after V _s	0 ms	

Pacing polarity

Parameter	Range of values	Standard
Polarity pace LV	LV tip/LV ring (bipolar 1) LV tip/RV ring (common ring bipolar 2) LV ring/LV tip (inverse bipolar 3) LV ring/RV ring (ring-to-ring bipolar 4)	LV tip/LV ring

Mode switching

Parameter	Range of values	Standard
Intervention rate	OFF 100 ... [10] ... 250 ppm	160 ppm
Activation criterion X	3 ... [1] ... 8	5
Deactivation criterion Y	3 ... [1] ... 8	5
Mode	VDD; VDD(R) when permanent VDD(R)	VDD
Ventricular pacing	RV; BiV	BiV

Parameter	Range of values	Standard
LV T-wave protection	OFF; ON	ON
Triggering	OFF RVs; RVs + PVC	RVs
Modification of basic rate	OFF; +5 ... [5] ... +30 ppm	+10 ppm

Automatic threshold monitoring

Automatic threshold monitoring (ATM)

Parameter	Range of values	Standard
ATM RV	ON; OFF	ON
ATM LV	ON; OFF	ON

Rate Adaptation

Accelerometer

Parameter	Range of values	Standard
Maximum sensor rate	AUTO 90 ... [5] ... 160 ppm	120 ppm
Sensor gain	1.0; 1.1; 1.3; 1.4; 1.6; 1.8; 2.0; 2.2; 2.6; 3.0; 3.3; 3.7; 4.0; 4.5; 5.0; 6.0; 7.0; 8.0; 9.0; 10; 11; 12; 14; 16; 18; 20; 22; 24; 28; 32; 35; 40	6.0
Automatic gain	OFF; ON	OFF
Sensor threshold	<ul style="list-style-type: none"> • Very low = 0 • Low = 3 • Medium = 7 • High = 11 • Very high = 15 	Medium
Rate increase	0.5; 1 ... [1] ... 6 ppm/cycle	2 ppm/cycle
Rate decrease	0.25 ... [0.25] ... 1.25 ppm/cycle	0.5 ppm/cycle

Tachyarrhythmia Therapy - Parameters

Sensing Parameters

Atrium

Lumax VR-T DX, DR-T; HF-T

Parameter	Range	Standard
Sensing	<ul style="list-style-type: none"> • STD - standard • OFF - disabled • IND - individual change in sensing details 	STD
Minimum threshold	0.2 ... [0.1] ... 2.0 mV	0,4 mV
Far-field protection after Vp	50 ... [25] ... 225 ms	75 ms
Far-field protection after Vs	OFF 25 ... [25] ... 225 ms	75 ms
Upper threshold DR, HF	50; 75; 87,5%	50%
Upper threshold DX	50; 75; 87,5%	75%
Lower threshold	12,5; 25; 50%	25%

Right ventricle

Lumax VR-T, VR-T DX, DR-T, HF-T

Parameter	Range	Standard
Sensing RV	<ul style="list-style-type: none"> • STD – standard • TWS – extended T-wave suppression • VFS – extended VF sensitivity • IND – individual change in sensing details 	STD
Minimum threshold	0.5 ... [0.1] ... 2.5 mV	0.8 mV
Blanking after atrial pacing	50 ... [10] ... 100 ms	50 ms
Upper threshold	50; 75; 87,5%	50%
Hold of upper threshold	100 ... [20] ... 600 ms	360 ms
Lower threshold	12,5; 25; 50%	25%

Parameter	Range	Standard
Blanking after RV pacing	AUTO 100 ... [10] ... 350 ms	AUTO
Post pace T-wave suppression	ON; OFF	OFF

Left ventricle

Lumax HF-T

Parameter	Range	Standard
Sensing LV	<ul style="list-style-type: none"> STD - standard OFF - disabled IND - individual change in sensing details 	STD
Minimum threshold	0.5 ... [0.1] ... 2.5 mV	1,6 mV
Upper threshold	50; 75; 87.5%	50%
Hold of upper threshold	100 ... [20] ... 600 ms	360 ms
Lower threshold	12,5; 25; 50%	50%

Polarity sense

Lumax HF-T

Parameter	Range	Standard
LV polarity sense	UNIP (LV tip/housing) BIPL (LV tip/LV ring)	UNIP

Detection Parameters

Detection intervals

Parameter	Range	Standard
Interval VT1	OFF 270 ... [10] ... 600 ms	OFF
Interval VT2	OFF 270 ... [10] ... 500 ms	OFF
Interval VF	OFF 200 ... [10] ... 400 ms	300 ms

Detection counter

Parameter	Range	Standard
Detection counter VT1	10 ... [2] ... 60	26
Detection counter VT2	10 ... [2] ... 40	16
Detection counter VF - X	6 ... [1] ... 30	8
Detection counter VF - Y	8 ... [1] ... 31	12

Onset

Parameter	Range	Standard
Onset criterion in VT1/VT2 with SMART	20%	20%
Onset criterion in VT1/VT2 without SMART	OFF 4 ... [4] ... 32%	20%

Stability

Parameter	Range	Standard
Stability in VT1/VT2 with SMART	12%	12%
Stability in VT1/VT2 without SMART	OFF 8 ... [4] ... 48 ms	24 ms
Stability ATP One Shot	12%	12%

SMART detection

Parameter	Range	Standard
SMART detection VT1	OFF; ON	ON
SMART detection VT2	OFF; ON	ON

Sustained VT; without SMART, without SMART redetection

Parameter	Range	Standard
Sustained VT	OFF 00:30; 01:00; 02:00; 03:00; 05:00; 10:00; 15:00; 20:00; 25:00; 30:00 [mm:ss]	OFF

Forced termination; with SMART, SMART redetection

Parameter	Range	Standard
Forced termination	OFF 1 ... [1] ... 15 min	1 min

Redetection counter

Parameter	Range	Standard
Redetection counter VT1	10 ... [2] ... 30	20
Redetection counter VT2	10 ... [2] ... 30	14

Therapy Parameters

Ventricle

Parameter	Range	Standard
Energy of 1st shock and 2nd shock VT1, VT2	OFF 1 ... [1] ... 16 ... [2] ... 40 J	40 J
Energy of 1st shock and 2nd shock VF	1 ... [1] ... 16 ... [2] ... 40 J	40 J
Shock sequences 3rd-nth shock in VT1	OFF 1*40 J; 2*40 J; 3*40 J; 4*40 J; 5*40 J; 6*40 J	6*40 J
Shock sequences 3rd-nth shock in VF	OFF 4*40 J; 5*40 J; 6*40 J	6*40 J
Number of shocks (VT1/VT2)	0 ... [1] ... 8	8
Number of shocks (VF)	2; 6 ... [1] ... 8	8
Confirmation in VT1, VT2, VF	OFF; ON	ON
Shock form in VT1, VT2, VF	Biphasic; biphasic2	Biphasic
Polarity in VT1, VT2, VF	Normal; inverse; alternating	Normal
Shock path (valid for all shocks including the painless shock impedance)	RV > SVC + housing RV > housing RV > SVC	RV > SVC + housing
Shock path DX		RV > housing

ATP parameters

Parameter	Range	Standard
ATP type	Burst; Ramp; Burst + PES	Burst
ATP attempts	OFF 1 ... [1] ... 10	3
Ventricular pacing (HF devices only)	RV; LV; BiV	RV
Number S1	1 ... [1] ... 10	5
Add. S1	OFF; ON	ON
R-S1 interval	200 ... [10] ... 500 ms (absolute); 70 ... [5] ... 95% (adaptive)	80%
S1 decrement	5 ... [5] ... 40 ms	10 ms
Scan decrement	OFF 5 ... [5] ... 40 ms	OFF
S1-S2 interval	200 ... [10] ... 500 ms (absolute); 70 ... [5] ... 95% (adaptive)	70%
Minimum ATP interval	200 ... [5] ... 300 ms	200 ms
ATP timeout	OFF 00:15 ... [00:15] ... 05:00 mm:ss	OFF
ATP optimization	OFF; ON	OFF
ATP pulse amplitude	7.5 V	7.5 V
ATP pulse width	1.5 ms	1.5 ms

ATP One Shot parameters; ATP in VF

Parameter	Range	Standard
ATP type	OFF; Burst; Ramp; Burst + PES	OFF
ATP attempts	1	1
Ventricular pacing (HF devices only)	RV; LV; BiV	RV
Number S1	1 ... [1] ... 10	8

Parameter	Range	Standard
R-S1 interval	200 ... [10] ... 350 ms (absolute); 70 ... [5] ... 95% (adaptive)	85%
S1 decrement	5 ... [5] ... 40 ms	10 ms
S1-S2 interval	200 ... [10] ... 350 ms (absolute); 70 ... [5] ... 95% (adaptive)	70%
Stability	12%	12%
ATP pulse amplitude	7.5 V	7.5 V
ATP pulse width	1.5 ms	1.5 ms

Progressive course of therapy

Parameter	Range	Standard
Progressive course of therapy	OFF; ON	ON

Post-shock pacing

Parameter	Range	Standard
Pacing mode (automatically programmed)	<ul style="list-style-type: none"> • DDI when permanent DDD(R), DDI(R), AAI(R) • VDI when permanent VDD(R), VDI(R) • VVI when permanent VVI(R); • OFF 	
Basic rate	30 ... [5] ... 100 ... [10] ... 160 ppm	60 ppm
Rate hysteresis	OFF -5 ... [-5] ... -65 ppm	OFF
AV delay	50 ... [10] ... 100 ... [10] ... 160 ms	140 ms
Ventricular pacing	RV; BiV	RV
LV T-wave protection	OFF; ON	ON
Triggering	OFF RVs; RVs + PVC.	OFF
Post-shock duration	OFF 00:10 ... [00:10] ... 00:50; 01:00 ... (01:00) ... 10:00 mm:ss	00:10 mm:ss

Diagnostics Functions

Home Monitoring, Holter Memory and Statistics

Home Monitoring

Parameter	Range	Standard
Home Monitoring	OFF; ON	OFF
Time of transmission	00:00 ... [1 min] ... 23:59 hh:mm	01:00 hh:mm
IEGM for therapy episodes	OFF; ON	OFF
IEGM for monitoring episodes	OFF; ON	OFF
Periodic IEGM	OFF 2; 3; 4; 6 months	OFF
Sustained atrial episodes (DR, HF)	OFF 0.5; 6; 12; 18 h	12 h
Thoracic impedance	ON; OFF	OFF

Holter episodes

Parameter	Range	Standard
AT/AF	OFF; ON	ON
SVT	OFF; ON	ON

Statistics

Parameter	Range	Standard
AT/AF rate	100 ... [10] ... 250 bpm	200 bpm
Start resting period	0:00 ... [10] ... 23:50 hh:mm	02:00 hh:mm
Resting period duration	0:00 ... [10] ... 23:50 hh:mm	04:00 hh:mm
Automatic impedance measurement	VR, DX, DR: ON; OFF HF: ON; OFF; LV OFF	ON
Thoracic impedance	ON; OFF	OFF

Characteristics

Electrical Characteristics

Measuring conditions

If not indicated otherwise, all specifications refer to the following conditions:

- Ambient temperature of $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ and $500\Omega \pm 5\%$ load (pacing/sensing) or $50\Omega \pm 1\%$ load (shock)
- Software used: Programmer in standard program

Standards

The specifications are made according to EN 45502-2-2:2008.

Factory settings

All therapy parameters are deactivated at delivery.

Arrhythmia classes VT1, VT2, VF	Antibradycardia pacing	Home Monitoring
OFF	OFF	OFF


Home Monitoring

Telemetry data

Nominal carrier frequency	Maximum power of transmission
403.62 MHz	< 25 μW (-16 dBm)


Atrial pacing pulse

Lumax DR-T; Lumax HF-T


Ventricular pacing pulse

Lumax VR-T; Lumax VR-T-DX; Lumax DR-T right ventricle; Lumax HF-T right and left ventricle; except for polarity pace LV ring/LV tip inverse bipolar

Lumax HF-T left ventricle; polarity pace LV ring/LV tip inverse bipolar

Resistance to interference

With unipolar sensing, the requirements of EN 45502-2-2:2008 are met for interference voltages < 0.6 mV (peak to peak).

EMC requirements relating to section 27.5 are met as long as atrial sensitivity is set to 1.0 mV (corresponds to factory settings) or higher (≥ 1.0 mV). Appropriate measures must be taken to assure interference-free therapy if more sensitive values are chosen.

Common mode rejection

Indication of rate in Hz, of common mode rejection in dB:

Rate	Atrium	Atrium	V right	V left
	VR-T DX	VR-T, DR-T, HF-T		
CMRR (Common Mode Rejection Ratio)				
16.6	67	65	64	64
50	80	80	67	64
60	81	75	66	63

ATP amplitude

Set parameters:

Burst	Value
Amplitude	7.5 V
Pulse width	1.5 ms
R-S1 interval	300 ms
Number S1	5
Load resistance	500 Ω

Results:

ATP amplitude	Value with tolerance	Measured minimum	Measured maximum
RV	7.5 V ± 1.5 V	7.42 V	7.50 V
Mean value	—	5.18 V	5.18 V
LV	7.5 V ± 1.5 V	7.50 V	7.50 V
Mean value	—	5.18 V	5.18 V
BiV	7.5 V ± 1.5 V	7.42 V	7.50 V
Mean value	—	5.16 V	5.18 V

Sensitivity setting

Indication of minimum and maximum of automatic sensitivity setting; measurement of the actual values with positive and negative polarity:

Note: The programmed atrial sensitivity is intensified by a factor of 5 for Lumax 540 VR-T DX.

Sensitivity	Test signal wave shape	Value	Tolerance	Measured value
A: positive [mV]	\sin^2 15 ms	0.2	0.2 ... 0.5	0.26
	Standard triangle	—	—	0.41
A: negative [mV]	\sin^2 15 ms	0.2	0.2 ... 0.5	0.28
	Standard triangle	—	—	0.43
VR-T DX: A: positive [mV]	\sin^2 15 ms	0.2	0.02 ... 0.10 (0.1 - 0.5)	0.08
	Standard triangle	—	—	0.08
VR-T DX: A: negative [mV]	\sin^2 15 ms	0.2	0.02 ... 0.10 (0.1 - 0.5)	0.10
	Standard triangle	—	—	0.12
RV: positive [mV]	\sin^2 40 ms	0.5	0.3 ... 0.7	0.56
	Standard triangle	—	—	0.58
RV: negative [mV]	\sin^2 40 ms	0.5	0.3 ... 0.7	0.45
	Standard triangle	—	—	0.58
LV: positive [mV]	\sin^2 40 ms	0.5	0.3 ... 0.7	0.64
	Standard triangle	—	—	0.65
LV: negative [mV]	\sin^2 40 ms	0.5	0.3 ... 0.7	0.54
	Standard triangle	—	—	0.68

Shock energies / peak voltages at: RV > SVC + housing

Lumax 540 VR-T, VR-T DX, DR-T, HF-T as 40 J device with shock path:
RV > SVC + housing

Shock energy set		Measurement of shock energy at 50 Ω
1 J	E nominal	94% ±25% 0.7 ... 1.18 J
	E actual	0.81 J
	Tolerance range: U max	90 ... 120 V
	U max actual	101.01 V
20 J	E nominal	94% ±10% 16.9 ... 20.9 J
	E actual	17.31 J
	Tolerance range: U max	440 ... 480 V
	U max actual	465 V
40 J	E nominal	94% ±10% 33.8 ... 41.4 J
	E actual	35.31 J
	Tolerance range: U max	630 ... 670 V
	U max actual	656.7 V

Shock energies / peak voltages at: RV > SVC

Lumax 540 VR-T, VR-T DX, DR-T, HF-T as 40 J device with shock path: RV > SVC

Shock energy set		Measurement of shock energy at 50 Ω
1 J	E nominal	94% ±25% 0.7 ... 1.18 J
	E actual	0.83 J
	Tolerance range: U max	90 ... 120 V
	U max actual	104.2 V

Shock energy set		Measurement of shock energy at 50 Ω
20 J	E nominal	94% ±10% 33.8 ... 41.4 J
	E actual	35.60 J
	Tolerance range: U max	630 ... 670 V
	U max actual	659.3 V
40 J	E nominal	94% ±10% 33.8 ... 41.4 J
	E actual	35.60 J
	Tolerance range: U max	630 ... 670 V
	U max actual	659.3 V

Behavior in case of EMI

In case of electromagnetic interference, the device switches into the VOO mode when the interference rate is exceeded if VVI[R] or VDD[R] pacing mode was set previously, or into the D00 pacing mode if DDI[R] or DDD[R] pacing mode was set previously.

Mechanical Characteristics

Housing

Device type	W x H x D in mm	Volume in ccm	Mass g
VR-T, VR-T DX, DR-T	66 x 55 x 13	37.2	92
HF-T	66 x 59 x 13	39.8	94

Materials in contact with body tissue

The following materials of the ICD system have contact with body tissue:

Housing	Header	Sealing plugs, sealing lips
Titanium	Epoxy resin	Silopren

Radiopaque marker

The radiopaque marker is on the left, positioned vertically.

Manufacturing year	Radiopaque marker
nn	BIO SH

ICD leads

Connections

Pacing/sensing	Shock
• VR-T, VR-T DX: 1 x IS-1 bipolar	• All device types: 2 x DF-1
• DR-T: 2 x IS-1 bipolar	
• HF-T: 3 x IS-1 bipolar	

ICD lead configuration

Connection	Lead connector
RV	DF-1 connector of the ventricular shock coil
SVC	DF-1 connector of the vena cava electrode
SA	IS-1 connector of the atrial sense electrode, bipolar
P/S A	IS-1 connector of the atrial pace/sense electrode, bipolar
P/S V	IS-1 connector of the ventricular pace/sense electrode, bipolar
CS	IS-1 connector of the ventricular pace/sense coronary sinus electrode, bipolar or unipolar

Tolerances

The following tolerances apply to the programmable parameters:

Parameter	Tolerance
Pacing pulse amplitudes	+20% to -25%
Pacing pulse width	±10%
Atrium sensitivity DR-T, HF-T	For values < 0,4 mV: 0.2 mV to 0.52 mV For values ≥ 0,4 mV: +30% to -50%
Atrium sensitivity VR-T DX	For values < 0,4 mV: -0.1 mV to +0.3 mV; for values > 0,4 mV: +30% to -60%
Sensitivity of right and left ventricle	±40%
Pacing intervals	±20 ms
Escape interval	±20 ms
Detection intervals	±20 ms
Refractory periods	±20 ms
ATP coupling interval	±20 ms

Parameter	Tolerance
Shock coupling interval	±20 ms
Shock energy/charging voltage	Electrical characteristics: see peak voltages

Included Equipment Lumax

Standard

In the sterile container there is:

- Device of the Lumax product family
- Blind plug DF, mounted in the header
- Blind plug IS-1, mounted in the header (only for HF devices)
- Screwdriver

In the storage package there is:

- Patient's ICD manual
- Sterile cover for the PGH programming head of the programmer
- Implantation protocol
- Serial number label
- Patient ID card
- Follow-up report
- Warranty card
- Technical manual

Orders

The devices are available with the following order numbers (some device types are not available in all countries):

Device type	Order number
540 VR-T	360348
540 VR-T DX	368852
540 DR-T	360346
540 HF-T	360347

Accessories

All accessories must correspond to the requirements of the EC Directive 90/385/EEC.

- BIOTRONIK leads
- BIOTRONIK programming and monitoring devices
- Permanent magnet M-50
- For Home Monitoring: BIOTRONIK transmitters

Service Life

Lumax Battery Data

Battery LiS R7

LITRONIK LiS 3192 R7:

Model	LiS 3192 R7
Manufacturer	LITRONIK GmbH & Co, 01796 Pirna, Germany
System	LiMnO ₂
Usable capacity until EOS	1720 mAh
Battery voltage at BOL (BOS)	3.2 V
Battery ID number shown on the programmer	3

Battery GB

GREATBATCH GB 2491:

Model	GB 2491
Manufacturer	GREATBATCH, INC. Clarence, NY 14031
System	Li/SVO/CFx
Usable capacity until EOS	1720 mAh
Battery voltage at BOL (BOS)	3.2 V
Battery ID number shown on the programmer	1

Lumax Battery Charge Times

Lumax 540 VR-T

Shock energy	Battery type	Order number	Charge time at BOL (BOS)	Charge time at ERI
HE	LiS 3192 R7 or GB 2491	360348	≤ 10 s	≤ 12 s

Lumax 540 VR-T DX

Shock energy	Battery type	Order number	Charge time at BOL (BOS)	Charge time at ERI
HE	LiS 3192 R7 or GB 2491	368852	≤ 10 s	≤ 12 s

Lumax 540 DR-T

Shock energy	Battery type	Order number	Charge time at BOL (BOS)	Charge time at ERI
HE	LiS 3192 R7 or GB 2491	360346	≤ 10 s	≤ 12 s

Lumax 540 HF-T

Shock energy	Battery type	Order number	Charge time at BOL (BOS)	Charge time at ERI
HE	LiS 3192 R7 or GB 2491	360347	≤ 10 s	≤ 12 s

Lumax Service Times

Legend

Legend for the following service time graphics:

- X axis: Number of maximum energy charging processes (shocks and capacitor reforming)
- Y axis: Longevity in years;
of Home Monitoring devices:
 - 1 periodic message per day
 - 12 IEGM transmissions per year
- Maximum energy charging processes per year:
 - A = 4
 - B = 8
 - C = 12
 - D = 16
 - E = 20

Calculation of service times


- The service times are calculated with the following values:
 - Pulse amplitude: 2.8 V
 - Pulse width: 0.4 ms
 - Pacing impedance: 500Ω
 - Basic rate: 60 ppm
- Pacing site:
 - Triple-chamber devices: RV
 - Dual-chamber devices: RV and A
 - Triple-chamber devices: RV, LV and A

Note: Devices in the Lumaxfamily should be implanted within 16 months between the date of manufacture and the use by date as specified on the package.

- If the ICD is implanted shortly before the use by date, the expected service time may be reduced by up to 15 months.
- Please note that four capacitor reformings are carried out per year. Therefore you have to account for at least 4 shocks per year, even if less than 4 have been delivered (see the following charts of service times).


Lumax 540 VR-T

Service times with LiS 3192 R7 and GB 2491 battery


Lumax 540 VR-T DX

Service times with LiS 3192 R7 and GB 2491 battery


Lumax 540 DR-T

Service times with LiS 3192 R7 and GB 2491 battery

Lumax 540 HF-T

Service times with LiS 3192 R7 and GB 2491 battery

Legend for the Label

The label icons symbolize the following:

	Manufacturing date
	Use by
	Storage temperature
	BIOTRONIK order number
	Serial number
	Product identification number
	Warning: dangerous voltages!
	Sterilized with ethylene oxide
	Do not resterilize!
	Single use only. Do not reuse!
	Non-sterile
	Follow the instructions in the technical manual!
	Contents
	Do not use if package is damaged!

	CE mark
	Transmitter with non-ionizing radiation on specified frequency; turned off during shipment
	Device with NGB coding and name of compatible ICD leads (example)
	Device with factory settings: Bradycardia and tachycardia therapies are set to OFF and/or the program settings are physiologically ineffective
	Screwdriver
	Position of connector ports in the header (example)
	Bipolar IS-1 connector
	Unipolar IS-1 connector
	Unipolar DF-1 connector
A	Atrium
V	Ventricle
CS	Coronary sinus
Pace	Pacing
Sense	Sensing
Shock	Shock
RV	High-voltage connector port
SVC	High-voltage connector port

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All rights reserved. Specifications subject to modification, revision and improvement.
® BIOTRONIK Home Monitoring, IEGM-Online HD and SMART Detection are registered trademarks of BIOTRONIK SE & Co. KG