# ScoutPro

## Kardiale Resynchronisationstherapie Cardiac Resynchronization Therapy

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Slitter Tool Advanced

## Table of Contents

<b>ScoutPro CS Lead Introducer System .....</b>	<b>22</b>
General Considerations .....	22
Indications .....	22
Contraindications .....	23
Potential Negative Side Effects .....	23
Scope of Delivery .....	24
Packaging and Sterility .....	26
Basic Safety Warnings .....	26
Handling Instructions .....	27
Disclaimer .....	32
Product Line .....	32
<b>Notes on Additionally Available Guiding Catheters .....</b>	<b>35</b>
Handling Instructions .....	37
<b>Legend for the Label .....</b>	<b>38</b>

# ScoutPro CS Lead Introducer System

## General Considerations

**Description** The BIOTRONIK ScoutPro CS (Coronary Sinus) lead introducer system is a set of complementary tools designed to be used together during the implantation of suitable leads in the coronary sinus.

Components of this lead introducer system as well as additional guiding catheters can also be ordered separately (see "Product Line", page 32 et seq.).

The ScoutPro lead introducer system is available in two nominal diameters: 7 French (ScoutPro 7 F) and 8 French (ScoutPro 8 F).

**Intended use** ScoutPro is suitable for use with approved polyurethane or polyurethane-coated leads with an appropriate diameter.

Lead introducer system	Suitable for leads with a diameter of less than:
ScoutPro 7 F	6.6 F
ScoutPro 8 F	7.8 F

Using ScoutPro has the following advantages:

- Quick probing of the ostium of the coronary sinus
- Easy access to the coronary sinus
- Monitoring of the coronary vessels using angiography
- Easy positioning and repositioning of the leads

## Indications

The ScoutPro CS lead introducer system and its separately available accessories are used to facilitate lead implantation in the left side of the heart via the coronary sinus.

## Contraindications

Use of the ScoutPro CS lead introducer system is contraindicated for the following:

- Patients with an existing or possible occlusion of the coronary vessels or inadequate anatomy of the coronary veins
- Patients with active systemic infection

## Potential Negative Side Effects

Possible complications are as follows:

- Allergic reactions to contrast media
- Hematoma
- Bleeding
- Infection
- Air embolism
- Pneumothorax
- Cardiac tamponade
- Myocardial damage
- Venous or cardiac perforation
- Local tissue response, fibrotic tissue formation
- Thromboembolism
- Damage to the heart valves
- Vascular occlusion
- Chronic nerve damage

## Scope of Delivery

### CS lead introducer system

**ScoutPro 7 F, order number: 350 239 or**

**ScoutPro 8 F, order number: 342 262**


Fig. 1: Package contents of the ScoutPro CS lead introducer system

## Individual items and technical data

Item	Technical data	ScoutPro 7 F	ScoutPro 8 F	
Puncture accessories				
1 Puncturing cannula		18 Gauge	18 Gauge	
1 Syringe	Volume:	10 ml	10 ml	
1 Guide wire	Length: Material:	Approx. 100 cm Stainless steel	Approx. 100 cm Stainless steel	
1 Peel-away introducer sheath	Length: Inner diameter: Material:	Approx. 14.3 cm 3.3 mm (10 F) PE	Approx. 14.3 cm 3.7 mm (11 F) PE	
Introducer accessories				
Guiding catheter	Working length: Soft tip length:  Inner diameter: Outer diameter: Internal material: External material:	<b>Sheath "MPEP"</b> 45 cm 9 mm + 1.7 mm (white soft tip)  2.38 mm (7.1 F)  3.0 mm (9.0 F) PTFE PEBAX	<b>Sheath "BIO 1"</b> 45 cm 12 mm + 2 mm (white soft tip)  2.8 mm (8.5 F)  3.4 mm (10.4 F) PTFE PEBAX	
	X-ray recognition:	Barium sulfate in the entire catheter area, distal marker ring (1 mm wide)		
Guiding catheter	Working length: Soft tip length:  Inner diameter: Outer diameter: Internal material: External material:	<b>Sheath "BIO 2"</b> 45 cm 9 mm + 1.7 mm (white soft tip)  2.38 mm (7.1 F)  3.0 mm (9.0 F) PTFE PEBAX	<b>Sheath "BIO 2"</b> 45 cm 12 mm + 2 mm (white soft tip)  2.8 mm (8.5 F)  3.4 mm (10.4 F) PTFE PEBAX	
	X-ray recognition:	Barium sulfate in the entire catheter area, distal marker ring (1 mm wide)		
1 Dilator (for the guiding catheters)	Working length: Outer diameter: Material:	Approx. 53 cm 2.3 mm [7 F] PE	Approx. 53 cm 2.7 mm [8 F] PE	
1 Hemostatic valve		With side connection and three-way spigot		
2 Slitter tools		Suitable for leads with a diameter of:		
		Slitter tool, green: Slitter tool, white:	4.1 F – 5.4 F 5.5 F – 7.2 F	4.1 F – 5.4 F 5.5 F – 7.2 F

Table 1: Scope of delivery and technical data

## Packaging and Sterility

The ScoutPro CS lead introducer system is delivered in a sterile plastic packaging consisting of a plastic container and cover (see Fig. 1). This plastic packaging is sealed in two sterile bags, one within the other, and sterilized with ethylene oxide, so that the inner bag is sterile on the outside.

Outer box and labeling

The sterile bags are packaged in an outer container that carries a quality control seal and a product information label. The label lists the model type, technical data, serial number, expiration date, and information about sterility and storage of the package.

- 1 To ensure sterility, the container should be checked for integrity prior to opening.
- 2 If a breach of sterility is suspected (open or damaged container) return the ScoutPro CS lead introducer system to BIOTRONIK.

**Note:** The implantation accessories of the ScoutPro CS lead introducer system are intended for single use only.

## Basic Safety Warnings

Not for multiple usage

The implantation accessories of the ScoutPro CS lead introducer system are intended for single use only.

Avoiding thrombosis

To inhibit blood clotting, BIOTRONIK recommends intravenously administering 50-100 IU heparin per kg body weight if the patient is not otherwise anticoagulated. The heparin may be directly administered through the ScoutPro guiding catheter or a different venous access point. In this case, please flush with sufficient physiological NaCl solution to guarantee that the heparin completely enters the patient's circulation.

Damage to the lead

Avoid damaging the lead with sharp-edged instruments - especially with the blade of the slitter tool when removing the catheter.

Contrast media

Only use water soluble contrast media.

## Handling Instructions

### Opening the packaging

- |                       |  |
|-----------------------|--|
| In a non-sterile area | Remove the outer sterile bag from the outer container and open it.         |
| In a sterile area     | Remove and open the inner sterile bag and then open the plastic packaging. |

### Preparation

- 1 Prior to use, all components should be thoroughly rinsed and, if possible, vented.
- Note:** Use a physiological NaCl solution for rinsing.
- 2 Choose one of the included or separately available ScoutPro guiding catheters with regard to the patients's atrial anatomy, as well as diameter, length, and curve form.
- 3 Prior to use, thoroughly rinse, ventilate, and remove the protective caps from the three-way spigot and connect the hemostatic valve to the selected guiding catheter.
- 4 Set the valve so that hemostasis is maintained during the complete procedure.

**Note:** The hemostatic valve can be tightened and loosened using the rotation ring.

**Note:** The three-way spigot allows 4 positions for different flow directions for contrast media, irrigation fluid, etc. indicated by the arrows in the figure.


Fig. 2: Possible positions of the three-way spigot:  
Position 2: Access to the valve is closed

- 5 Insert the associated dilator into the corresponding guiding catheter. This causes the catheter to straighten.

**Note:** Do not attempt to reshape the guiding catheter! It is preformed and resumes its curved shape after the dilator has been removed.

### Puncture

- 1 Connect the body of the syringe to the puncturing cannula. Puncture the selected vein (e.g., subclavian vein, cephalic vein) at a suitable site.

**Caution!** For patients with chronic lung disease, there is an increased risk of pneumothorax!

- 2 To check the venous puncture, aspirate blood into the syringe.
- 3 Remove the syringe body. Leave the cannula in position.

### Guide wire

- 1 Using X-ray detection, introduce the guide wire through the cannula into the vein and advance it up to the atrium.

Depending on the patient's body size, approximately 60 cm of the wire usually remains outside of the body.

- 2 Remove the cannula.

In addition, you can use the peel-away sheath with dilator by:

- Pushing the introducer, with the dilator fully in place, via the guide wire into the vessel, and then
- remove the dilator.  
The introducer sheath remains in the vessel.

- 3 Insert the selected guiding catheter with dilator and mounted hemostatic valve using the guide wire. Advance the guiding catheter up to the atrium.

**Caution!** Carefully advance the guiding catheter to prevent damage to the vessel walls!

- 4 Remove the dilator and the guide wire. The guiding catheter resumes its preshaped form, facilitating access to the coronary sinus.
- 5 Advance the guiding catheter through the ostium and into the coronary sinus.

### Probing the ostium

Probing of the ostium can be simplified with the use of a suitable catheter.

The AL-1 and MPA 2 diagnostic catheters are suitable for the ScoutPro CS lead introducer system.

**Note:** Follow the instructions in the appropriate manual when using the diagnostic catheter AL-1 or MPA 2.

- 1 Open the valve using the rotation ring in order to insert the probing catheter.
- 2 Advance the tip of the respective probing catheter through the guiding catheter into the coronary sinus.
- 3 Advance the guiding catheter into the coronary sinus using the probing catheter.
- 4 Remove the probing catheter.

### Angiography

A coronary angiography can be performed to facilitate lead positioning.

- This is performed by using a venogram balloon catheter and inserting it into the coronary sinus using the guiding catheter.

**Note:** Please consult the user manual provided with the venogram balloon catheter.

### Inserting the lead

- Advance the lead tip through the guiding catheter until a suitable site is reached.

**Note:** The ScoutPro CS lead introducer system (including all accessories) may only be used with approved polyurethane or polyurethane-coated leads.

### Subsequent steps: Removing the guiding catheter with the slitter tool

**Note:** The guiding catheters from the BIOTRONIK ScoutPro CS lead introducer system cannot be lacerated by hand using the peel-away technique. These catheters must be slit open using a slitter tool.

Instead of the included slitter tools, you can also use the separately available Slitter Tool Advanced. Please then follow the included instructions for use of the accessory.

- 1 Dismantle the hemostatic valve and remove it over the proximal end of the lead.
- 2 Select the appropriate slitter tool from the ScoutPro CS lead introducer system or from other appropriate separately available slitter tools (see the last line in table 1 on page 25).
- 3 Attach the slitter tool to the lead near the proximal end of the guiding catheter.
- 4 For fixation, press the lead into the fixation groove of the slitter tool.


Fig. 3: Press the lead into the slitter tool

- 5 Tear off the peel-away proximal portion from the guiding catheter.
- 6 Hold the slitter tool and slide it along the guiding catheter, slitting open the guiding catheter with the steel blade of the slitter tool anchored to the lead.

- 7 If necessary, remove the peel-away introducer by hand or using the slitter tool.

**Note:** When using the slitter tool, make sure that the guiding catheter moves straight against the slitter tool (Fig. 4). Avoid cross threading the slitter tool since this may toughen the process. Moreover, the slitter tool could get stuck in the guiding catheter (Fig. 5). If, however, the slitter tool becomes stuck, replace it with a new slitter tool suitable for the lead you are using (same color; sold separately) or with the Slitter Tool Advanced (separately available) and continue slitting with the new tool.

Fig. 4: Correct handling: Guiding catheter is pulled straight against the slitter tool.

Fig. 5: Please avoid: An incorrectly guided slitter tool can get stuck in the catheter!

## Disposal

Used ScoutPro implantation accessories can be properly disposed of together with contaminated hospital refuse. They do not contain any materials that require special waste disposal.

## Disclaimer

BIOTRONIK's ScoutPro CS lead introducer system (including all individual accessories available separately) that is used in conjunction with devices for implantation in the coronary sinus have been qualified, manufactured, and tested according to proven and accepted standards and procedures. However, the physician must be aware that the accessories may easily be damaged through improper handling or use. Except as set forth in BIOTRONIK's limited warranty, BIOTRONIK makes no expressed or implied warranties for its accessories.

## Product Line

The ScoutPro line consists of the CS lead introducer system described in this technical manual and a number of additional, solely separately available accessories.

All components of the CS lead introducer system can be separately ordered and delivered, except for two.

The table gives a complete overview of our product line (see next page).

Item	Order number	in the Scout Pro system Separately available	
ScoutPro 7 F CS lead introducer system	350 239		
ScoutPro 8 F CS lead introducer system	342 262		
<b>Guiding catheter for ScoutPro 8 F</b>			
ScoutPro sheath BIO 1 (8 F)		x	
ScoutPro sheath BIO 2 (8 F)		x	
ScoutPro sheath Amplatz 6.0 (8 F)	342 264		x
ScoutPro sheath Multipurpose Hook (8 F)	342 265		x
ScoutPro sheath Hook (8 F)	342 263		x
<b>Guiding catheter for ScoutPro 7 F</b>			
ScoutPro 7 F sheath MPEP	359 373	x	x
ScoutPro 7 F sheath BIO 2	359 374	x	x
ScoutPro 7 F sheath Amplatz 6.0	350 235		x
ScoutPro 7 F sheath Multipurpose Hook	350 237		x
ScoutPro 7 F sheath Hook	350 236		x
ScoutPro 7 F sheath Extended Hook	359 369		x
ScoutPro 7 F sheath Extended Hook Right	359 370		x
ScoutPro 7 F sheath MPEP L	361 533		x
ScoutPro 7 F sheath BIO 2 L	361 529		x
ScoutPro 7 F sheath Amplatz 6.0 L	361 527		x
ScoutPro 7 F sheath Multipurpose Hook L	361 534		x
ScoutPro 7 F sheath Hook L	361 532		x
ScoutPro 7 F sheath Extended Hook L	361 530		x
ScoutPro 7 F sheath Extended Hook Right L	361 531		x
<b>Compatible accessory</b>			
Slitter tools "ScoutPro Slitter Tools"	345 963	x	x
Hemostatic valve "ScoutPro Hemostatic Valve"	345 968	x	x
Guide wire "ScoutPro Guidewire"; Stainless steel; approx. 100 cm long	359 805	x	x
Slitter tools "ScoutPro Slitter Tool Advanced"	356 552		x

Table 2: Product line of the ScoutPro 7 F and ScoutPro 8 F CS lead introducer systems

## Notes on Additionally Available Guiding Catheters

These ScoutPro guiding catheters can be separately ordered and delivered. Only the 7-F-catheter sheath "MPEP" and the sheath "BIO 2" are supplied in the CS lead introducer system ScoutPro.

- ScoutPro 7 F sheath "Hook"
- ScoutPro 7 F sheath "Multipurpose Hook"
- ScoutPro 7 F sheath "Amplatz 6.0"
- ScoutPro 7 F sheath "MPEP"
- ScoutPro 7 F sheath "BIO 2"
- ScoutPro 7 F sheath "Extended Hook"
- ScoutPro 7 F sheath "Extended Hook Right"
- ScoutPro 7 F sheath "Hook L"
- ScoutPro 7 F sheath "Multipurpose Hook L"
- ScoutPro 7 F sheath "Amplatz 6.0 L"
- ScoutPro 7 F sheath "MPEP L"
- ScoutPro 7 F sheath "BIO 2 L"
- ScoutPro 7 F sheath "Extended Hook L"
- ScoutPro 7 F sheath "Extended Hook Right L"
- ScoutPro 8 F sheath "Hook"
- ScoutPro 8 F sheath "Multipurpose Hook"
- ScoutPro 8 F sheath "Amplatz 6.0"

## Intended use

These guiding catheters may only be used in connection with the respective Scout Pro CS lead introducer system 7 F or ScoutPro 8 F or with the accessories ("ScoutPro Guidewire", "ScoutPro Hemostatic Valve", "ScoutPro Slitter Tool" and "ScoutPro Slitter Tool Advanced"), which can be ordered separately.

In addition to the two catheters included in the CS lead introducer system ScoutPro, there are further catheter lengths and curve forms available.

Select the appropriate catheter according to the patient's atrial anatomy.

The guiding catheters are suitable for use with approved polyurethane or polyurethane-coated leads with an appropriate diameter.

Guiding catheter for the CS lead introducer system	Suitable for leads with a diameter of less than:
ScoutPro 7 F	6.6 F
ScoutPro 8 F	7.8 F

## Indications, contraindications, potential complications, and basic safety warnings

The indications, contraindications, and potential complications (negative side effects) of the ScoutPro lead introducer system (see page 22 and 23) also apply to the additional guiding catheters. All basic safety warnings in this technical manual apply correspondingly to the use of the additional guiding catheters.

## Packaging and sterility

Each guiding catheter is sterile and supplied with a dilator. They are packaged in two sterile bags, one within the other, and sterilized with ethylene oxide, so that the inner bag is also sterile on the outside.

Outer box and  
labeling

The sterile bags are packaged in an outer container that carries a quality control seal and a product information label. The label lists the model type, technical data, serial number, expiration date, and information about sterility and storage of the package.

- 1 Examine the package for damage before opening to ensure sterility.
- 2 If a breach of sterility is suspected (open or damaged container), return the ScoutPro CS lead introducer system to BIOTRONIK.

**Note:** The guiding catheters are intended for single use only.

## Handling Instructions

### Opening the packaging

In a non-sterile area	Remove the outer sterile bag from the outer container and open it.
In a sterile area	Remove and open the inner sterile bag.

### Using the guiding catheter with the ScoutPro CS lead introducer system

The handling of the additional guiding catheters in connection with the ScoutPro CS lead introducer system or with the separately available accessories (ScoutPro Guide Wire, ScoutPro Slitter Tools, and ScoutPro Hemostatic Valve) is outlined in section "Handling Instructions", page 27.

### Technical data

ScoutPro 7 F			ScoutPro 8 F
Guiding catheter	Models "..."	Models "... L"	
Working length	45 cm	50 cm	45 cm
Inner diameter	2.38 mm (7.1 F)		2.8 mm (8.5 F)
Inner material:	PTFE		
Outer diameter:	3.0 mm (9 F)		3.4 mm (10.4 F)
External material:	PEBAX		
X-ray contrast	Barium sulfate in the entire catheter area		
X-ray marker on the distal end	Marker ring (width: 1 mm)	Marker zone (width: 2 mm)	Marker ring (width: 1 mm)
Soft tip length:	9 mm + 1.7 mm (white soft tip)		12 mm + 2 mm (white soft tip)
Dilator			
Working length:	53 cm	59 cm	53 cm
Material:	PE		
Outer diameter:	2.3 mm (6.9 F)		2.7 mm (8 F)

Table 3: Guiding catheter technical data

## Legend for the Label

Date of manufacture of the device

Expiration date; only valid if sterile container is undamaged!

Never use expired devices!

Acceptable temperature ranges for storage

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All rights reserved. Specifications subject  
to modification, revision and improvement.

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This product conforms with the EC  
directive 90/385/EEC relating to active  
implantable medical devices. It was  
approved by an independent Notified Body  
and is therefore designated with the CE  
mark. The product can be used in all  
European Union countries as well as in  
countries that recognize the above-  
mentioned directive.

 0123 (2003)  
90/385/EEC

 0123 (2005)  
90/385/EEC

 0123 (2007)  
90/385/EEC

 0123 (2009)  
90/385/EEC

09-D-44

Revision: H (2009-09-30)

350862

Manufacturer:  
BIOTRONIK SE & Co. KG  
Woermannkehe 1  
12359 Berlin · Germany  
Tel +49 (0) 30 68905-0  
Fax +49 (0) 30 6852804  
sales@biotronik.com  
www.biotronik.com

Sales organization in the USA:  
BIOTRONIK, Inc.  
6024 Jean Road  
Lake Oswego, OR 97 035-5369  
Tel (800) 547 -0394 (24-hour)  
Fax (503) 635 -9936  
marketing@biotronicusa.com

 **BIOTRONIK**