Philos II

DDDR, Dual-Chamber Pacemaker with BIOTRONIK Home Monitoring®
DDDR, Zweikammer-Schrittmacher mit BIOTRONIK Home Monitoring®
DDDR, marcapasos bicameral con BIOTRONIK Home Monitoring®
Pacemaker bicamerale DDDR con BIOTRONIK Home Monitoring®
植入式心脏起搏器 (DDDR, 带BIOTRONIK Home Monitoring® 的双腔起搏器)

Technical Manual • en
Gebrauchsanweisung • de
Manual técnico • es
Manuale tecnico di istruzione • it
技术手册 • zh

BIOTRONIK
excellence for life

en	• English	2
de	• Deutsch	20
es	• Español	38
it	• Italiano	56
zh	• 中文	74

Product Overview

Intended Medical Use

Philos II is a family of pacemakers that may be used for all indications of bradycardia arrhythmia. The product family consists of 6 pacemakers. The Home Monitoring function allows external monitoring.

Implantation of a permanent pacemaker is generally indicated as therapy for symptomatic bradycardia.

Clinical symptoms include:

- Insufficient cerebral perfusion accompanied by presyncope and syncope
- Acute bouts of dizziness
- Chronic symptoms of reduced cardiac performance
- Uncharacteristic exhaustion or problems with concentrating

Device System

The Philos II device system includes the following parts:

- Single-chamber or dual-chamber device with connections for unipolar or bipolar sensing and pacing
- Leads
- ICS 3000 programmer with current device program

Device

The device's housing is made of biocompatible titanium, welded from outside and thus hermetically sealed. It serves as an antipole for unipolar lead configuration. The ellipsoid shape facilitates ingrowth into the pectoral muscle area. The labeling provides information about the device type and arrangement of the connections.

Leads

The leads are sheathed with biocompatible silicone. They can be flexibly maneuvered, are long-term stable, and are equipped for active or passive fixation. They are implanted using a lead delivery set. Some leads are coated with polyurethane which is known to increase the sliding properties for the lead.

The coating of steroid-eluting leads reduces inflammatory processes. The fractal design of the leads provides for low pacing thresholds, high pacing impedance, and a low risk of oversensing.

Programmer

The transportable ICS 3000 Implant Control System contains the device programs. It acts as a programmer, ECG and IEGM monitor and recorder, miniclinic, and storage medium.

The programmer communicates wirelessly with the device via the PGH programming head. The operation module of the programmer has a TFT touch screen with color display, on which the ECG, IEGM, marker and functions are shown simultaneously.

The ICS 3000 provides the following functions:

- Measures the pacing impedance
- Displays and prints real-time and saved IEGMs with annotated markers
- Determines the pacing threshold as well as several non-invasive induction techniques

Device/lead connection scheme

- BIOTRONIK pacemakers have been designed for use with leads with a unipolar or bipolar IS-1 connection.
- The IS-1 pacemaker lead connection system is in compliance with ISO 5841-3: 2000-10.

Connection scheme of the device types:

SR	S	DR(-T)	D	SLR
WIR/AAIR	VVI/AAI	DDDR	DDD	VDDR
				
IS-1	IS-1	IS-1	IS-1	IS-1

Device and lead connections:

	S, SR	D, DR-T	SLR
Atrium	IS-1 unipolar or bipolar	IS-1 unipolar or bipolar	IS-1 bipolar
Ventricle		IS-1 unipolar or bipolar	IS-1 unipolar or bipolar

BIOTRONIK Home Monitoring®

In addition to effective pacemaker therapy, a complete therapy management system is available.

To enable early detection of asymptomatic events and timely intervention if necessary, a pacing system with Home Monitoring checks its status at specific intervals, and provides information on the cardiac condition of the patient.

The CardioMessenger is a mobile modem with a broadband power supply for international application. It has a lithium ion battery that lasts at least 24 hours. The device transmits diagnostic, therapeutic, and technical information to the mobile transmitter CardioMessenger via the antenna in the header and using telemetry.

The data is sent to the BIOTRONIK Service Center by means of the cellular network. There, they are accessible to the patient's physician via the Internet.

Technical manuals

The following technical manuals provide information about the application of the pacemaker system:

- Philos II family, technical manual for the device
- ICS 3000 SW, technical manual for the device program, on CD and as online help
- Pacemaker leads, technical manual for the leads
- ICS 3000, technical manual
- ICS 3000 with Implant Module, supplement to the technical manual

Product Models

The Philos II product family consists of six pacemakers.

Philos II DR and D

The Philos II DR and Philos II D dual-chamber systems are used with separate atrial and ventricular leads. They are suitable for AV sequential pacing.

Philos II DR-T

The Philos II DR-T dual-chamber system is technically identical to the Philos II DR. Additionally, it is equipped with the Home Monitoring function.

Philos II SLR

The Philos II SLR dual-chamber system only requires a single lead. This single-lead system facilitates atrial synchronized ventricular pacing.

Philos II SR and S

The Philos II SR and Philos II S single-chamber systems only require one lead. They are suitable for therapy in the atrium or the ventricle.

NBG Code

The following NBG codes are according to Bernstein et al., The Revised NASPE/BPEG Generic Code for Antibradyardia, Adaptive-Rate, and Multisite Pacing. PACE 2002, Vol. 25, No. 2: 260-264.

Philos II DR, DR-T

D	Pacing in both chambers
D	Sensing in both chambers
D	Pulse inhibition and triggering
R	Rate adaptation

Philos II D

D	Pacing in both chambers
D	Sensing in both chambers
D	Pulse inhibition and triggering

Philos II SLR

V	Ventricular pacing
D	Sensing in both chambers
D	Pulse inhibition and triggering
R	Rate adaptation

Philos II SR

A/V	Pacing in both chambers
A/V	Sensing in both chambers
I	Pulse inhibition in A/V
R	Rate adaptation

Philos II S

A/V	Pacing in both chambers
A/V	Sensing in both chambers
I	Pulse inhibition in A/V

Therapeutic and Diagnostic Functions

Overview

Depending on the type, the device program contains all pacemaker functions for 1 or 2 chambers. All the systems have extensive features that allow quick diagnosis and delivery of safe therapy for cases of bradycardia arrhythmia. The follow-up assistant functions have been automated except for the retrograde conduction test. This saves the physician time, and eliminates problems in verifying and adjusting the pacemaker.

ACC – Active Capture Control

- The amplitude control function (ACC – Active Capture Control) continuously monitors the effectiveness of ventricular pacing and continuously adjusts the pacing amplitude to the threshold.

Automatic threshold measurement

- The ventricular threshold test can be performed automatically, because the ACC function measures the threshold periodically.

Automatic lead check

- The regular automatic lead impedance check triggers a switch from bipolar to unipolar pacing mode when values outside the normal range occur.

Timing

- Innovative rate hystereses promote the patient's intrinsic rhythm and avoid unnecessary overdrive pacing.
- AV hysteresis functions support intrinsic conduction and thus the natural contraction sequence.
- The night program adjusts the pacing rate to the reduced metabolic needs of the resting patient during the night.
- The Rate Fading (rate smoothing) function ensures that the heart rate does not drop abruptly during a spontaneous intrinsic decrease, but rather is gradually reduced until the basic or sensor rate is reached.

Rate adaptation

- The rate is adapted by means of an integrated motion sensor (R mode).

Antitachycardia functions

- Antitachycardia functions significantly protect the patient from the consequences of tachycardias. Automatic mode switching prevents atrial-controlled ventricular pacing in the event of atrial tachycardias.
- Comprehensive algorithms assist in the prevention, detection and termination of pacemaker-mediated tachycardia (PMT).

IEGM recordings

- IEGM recordings provide insight into the events before a tachycardic phase.

Diagnostic functions (statistics)

- Comprehensive memory functions (histogram, rate trend, activity chart) facilitate the evaluation of both patient and pacemaker status.
- Intrinsic AV conduction statistics are suitable for optimizing the programmed AV delay and AV hysteresis.
- Atrial and ventricular extrasystoles as well as atrial tachycardias can be analyzed and classified with respect to their complexity and timing (time of occurrence).

NIPS

- An external pulse control function (NIPS: non-invasive programmed stimulation) is available to terminate atrial tachycardias and for electrophysiologic examinations. Burst stimulation with real-time control of the burst rate and programmed stimulation up to 4 extrastimuli are available.
Backup pacing in the ventricle is always available while using NIPS in case of loss of ventricle contraction.

Simplified follow-up

- Automatic functions and the storage of follow-up data in the device simplify and accelerate the follow-up process.

Home Monitoring: message retrieval

This function is available in '-T' devices (e.g. DR-T):

The device sends information to the Service Center once a day and additional messages after certain events. All data can be viewed on the Internet subsequently.

Three types of messages are sent:

- Trend message: the time (daily) triggers the message.
- Event message: the event triggers the message.
- Patient message: the pacemaker patient triggers the message using a special magnet.

Prior to Implantation

General Indications

The following conditions apply to pacemakers:

- Sinus node arrest and symptomatic bradycardia with or without AV conduction disorder
- Intermittent or complete AV block
- Brady/tachycardia syndrome or other symptoms of sick sinus syndrome that result in symptomatic bradycardia
- Supraventricular reentry tachycardias that can be suppressed by chronic AV sequential pacing
- Atrial and ventricular ectopic arrhythmias that can be suppressed by permanent AV sequential pacing

Special Indications

Indications: dual-chamber pacemakers

The following conditions apply as special indications for dual-chamber pacemakers:

- In contrast to a single-chamber pacemaker, a dual-chamber pacemaker is indicated for patients who require increased cardiac output. This includes active patients and patients who have experienced, or are likely to experience, pacemaker syndrome.
- Atrial-controlled, dual-chamber modes (DDD, VDD) are indicated for patients who have an intact spontaneous atrial rhythm. Ventricular-controlled, AV sequential dual-chamber pacing modes (DDI, DVI, VDI) are indicated for patients for whom ventricular pulse triggering by spontaneous atrial events is not required or desired.

Rate-adaptive pacing

Special indications for rate-adaptive pacing are the following:

- Patients who exhibit chronotropic incompetence and require increased pacing rates with physical activity.

Single-chamber pacemakers

Special indications for single-chamber pacemakers are the following:

- The AAI mode is indicated in the presence of symptomatic sinus node dysfunction as long as adequate AV conduction exists.
- The VVI mode is indicated in cases of symptomatic bradyarrhythmia as long as there is no (longer) significant atrial contribution to hemodynamics.

Indications for other modes

Special indications for other modes are the following:

- Other modes as well as their asynchronous DDO, AOO, and VOO modes (when the sensing function is reduced) are indicated in case of medical/technical complications (e.g., electromagnetic interference, sensing defects, lead fracture, detection of myopotentials, muscle stimulation, etc.).
- The triggered pacing modes DDT, DDI/T, VDT, DVT, AAT, and VVT as well as the VDI and OFF modes are indicated for diagnostic purposes.

General Contraindications

There are no known contraindications for the use of multiprogrammable and multi-functional dual-chamber pacemakers, provided that implantation is preceded by an adequate diagnosis, and no parameter combinations inappropriate for the patient's condition are programmed.

Note: In individual cases, it is recommended that the tolerance and effectiveness of parameter combinations are checked by observing the patient for some time after programming.

Special Contraindications

Operating modes with atrial control

The following are contraindicated:

- Operating modes with atrial control (DDD, VDD, AAI) are contraindicated with chronic atrial tachycardia as well as chronic atrial fibrillation or flutter.

Retrograde conduction

- If slow retrograde conduction is encountered after ventricular pacing, a longer atrial refractory period and/or a shorter AV delay may have to be programmed to prevent pacemaker-mediated tachycardia. Programming to DDI, DVI, or VVI mode is rarely required in these instances.

High pacing rates

- If elevated rates above the basic rate are not well tolerated by the patient (e.g., angina pectoris) a low upper rate and low maximum sensor rate should be programmed. In these cases, atrial-controlled modes and rate-adaptive modes may even be contraindicated.

Pacemaker syndrome

- If a pacemaker syndrome has been observed or is likely to develop, the modes VDD, VVI, and VOO are contraindicated. The DDI mode is contraindicated in cases of pacemaker syndrome where sinus rates are higher than the basic rate.

Atrial single-chamber pacing

- Atrial single-chamber pacing is contraindicated in the presence of existing AV conduction disorders or if failing AV conduction can be demonstrated.

Competitive stimulation

- In cases of dominating intrinsic rhythms, modes without sensing and inhibiting functions are contraindicated to prevent competitive stimulation.

Unipolar pacing

- Unipolar pacing is contraindicated for patients who also have an implanted cardioverter-defibrillator (ICD). There is a risk of ICD inhibition or accidental delivery of pacemaker pulses.

Shipping and Storage

Ambient conditions for shipping and storage

Adhere to the ambient conditions for shipping and storage.

CAUTION

The device should only be shipped and stored at temperatures between 5°C and 55°C (41°F to 131°F). Exposure to temperatures outside this range may result in device malfunction.

Storage package

The storage package consists of the following components:

- Folded cardboard box with a quality control seal and a label with the following information:
 - Name of model
 - Technical details
 - Serial number
 - Use by date
 - Details on sterility
 - Information on storing the device
- The folded cardboard box contains the plastic container with the device and the documentation material.

Implantation

Sterility

Delivery

The device and the accessories have been gas sterilized. Sterility is guaranteed only if the plastic container and quality control seal have not been damaged.

Sterile container

The device and accessories are packaged respectively in two separately sealed plastic containers. The inner plastic container is also sterile on the outside so that it can be transferred in a sterile state during implantation.

Unpacking the device

Proceed as follows:

Step	Action
1	Peel off the sealing paper of the non-sterile outer plastic container at the marked position in the direction indicated by the arrow. The inner plastic container may not come into contact with persons who have not sterilized their hands or gloves, nor with non-sterile instruments.
2	Take hold of the inner plastic container by the gripping tab and take it out the outer plastic container.
3	Peel off the sealing paper of the sterile inner plastic container at the marked position in the direction indicated by the arrow.

Connecting the Leads

Precautionary measures

Note: Use only adapters approved by BIOTRONIK for leads with different connections. Should you have questions about the compatibility of other manufacturers' leads, please contact BIOTRONIK.

WARNING

Short circuit due to open lead connections

Open and thus not electrolyte-tight IS-1 connections may cause undesired current flow to the body and penetration of bodily fluid into the device.

- Close IS-1 connections that are not in use with IS-1 blind plugs.

CAUTION

When connecting unipolar leads to the pacemaker, you must set the pacing and/or sensing function of the respective channel to unipolar configuration.

Connecting the lead connector to the device

Leads are connected to the device according to the diagram on the header. Proceed as follows:

Step	Action
1	Disconnect stylets and insertion aids from the lead connector.
2	<ul style="list-style-type: none">Connect the unipolar or bipolar IS-1 lead connector atrium to A.Connect the unipolar or bipolar IS-1-lead connector ventricle to V.
3	Push the lead connector into the header without bending the conductor until the connector tip becomes visible behind the set screw block.
4	If the lead connector cannot be inserted completely, the set screw may be protruding into the cavity of the set screw block. Carefully loosen the set screw without completely unscrewing it, so that it does not become tilted upon retightening.
5	Use the screwdriver to perpendicularly pierce through the slitting in the center of the silicone plug until it reaches the set screw.
6	Turn the set screw clockwise until the torque control starts (you will hear a clicking sound).
7	Carefully withdraw the screwdriver without retracting the set screw. <ul style="list-style-type: none">When you withdraw the screwdriver, the silicone plug automatically seals the lead connection safely.

Implanting

Implantation site

In general the pacemaker is implanted subcutaneously or subpectorally on the right depending on the lead configuration as well as the anatomy of the patient.

Device status prior to implantation

The device is shipped with factory settings and can be connected and implanted in this state.

Sequence

Proceed as follows:

Step	Action
1	Shape the device pocket and prepare the vein.
2	Implant the leads and perform measurements.
3	Connect the device and leads.
4	Insert the device.
5	Guide the fixation suture through the opening in the header and fixate the device in the prepared pocket.
6	Close the device pocket.

Note: During auto-initialization (30 min.), the device cannot be re-programmed nor can a follow-up be conducted; otherwise, auto-initialization will be deactivated and will no longer be available.

Note: Perform a follow-up to check whether the device is functioning properly.

WARNING

Inadequate parameters or interrupted telemetry

The patient could enter a hemodynamically critical state due to the application of inadequate parameters or interruptions in telemetry while a temporary program controlled by the programmer is running.

- Continuously monitor the ECG and the patient's condition.
- Lift the programming head at least 30 cm. This serves to reactivate the permanent program.

After Implantation

Follow-up

Follow-up intervals

Follow-ups must be performed at regular, agreed intervals.

- Following the lead ingrowth phase, approximately 3 months after implantation, the first follow-up should be carried out by the physician using the programmer (in-office follow-up).
- The next in-office follow-up should be carried out once a year and no later than 12 months after the last in-office follow-up.

Follow-up with BIOTRONIK Home Monitoring®

Monitoring using the Home Monitoring function does not serve to replace regular in-office appointments with the physician required for other medical reasons. Follow-up supported by Home Monitoring can be used to functionally replace in-office follow-up under the following conditions:

- The patient was informed that the physician must be contacted if symptoms worsen or if new symptoms arise despite the use of the Home Monitoring function.
- Device messages are transmitted regularly.
- The physician decides whether the data transmitted via Home Monitoring with regard to the patient's clinical condition as well as the technical state of the device system are sufficient. If not, an in-office follow-up has to be carried out.

Possible early detection due to information gained via Home Monitoring may necessitate an additional in-office follow-up. For example, the data may indicate an early stage lead problems or a foreseeable end of service time (ERI). Furthermore the data could provide indications of previously unrecognized arrhythmias or modification of the therapy by reprogramming the device.

Follow-up with the programmer

Use the following procedure for in-house follow-up:

1	Record and evaluate the external ECG.
2	Check the sensing and pacing function.
3	Interrogate the device.
4	Evaluate the status and automatically measured follow-up data.
5	Possibly evaluate statistics and Holter/IEGM recording.
6	Manually perform standard tests if necessary.
7	Possibly customize program functions and parameters.
8	Transmit the program permanently to the device.
9	Print and document follow-up data (print report).
10	Finish the follow-up for this patient.

Elective Replacement Indication - ERI

Defined operating statuses of the device:

BOS	Beginning of Service	Battery is in good condition; normal follow-up.
ERI	Elective Replacement Indication	The replacement time has been reached. The pacemaker must be replaced.
EOS	End of Service	End of service time with regular pacemaker activity.

Determinants

The length of the period from beginning of service (BOS) until elective replacement indication (ERI) depends on several factors. This includes:

- Battery capacity
- Lead impedance
- Pacing program
- Pacing to inhibition ratio
- Pacemaker circuit properties

Activation

ERI detection is automatically activated if one of the following conditions occurs:

- Successfully executed auto-initialization.
- The programmer measures a lead impedance smaller than 3200 Ohm.
- The pacemaker is stored for longer than 24 months.
- The safe program is transmitted.

ERI display

- The pacemaker signals the replacement indication by means of a defined decrease both in the programmed basic rate and the magnet rate.
- The elective replacement indication ERI is displayed by the programmer upon interrogating the pacemaker.
 - In the course of the follow-up, the pacemaker displays the expected value until reaching ERI, based on the permanent program. This value is derived from the measured energy consumption of the battery.
 - A respective note with the obtained data can be printed out. The battery status can also be tested using analog telemetry.
 - If program parameters are modified, the calculated time until ERI for the edited program is also displayed on the programmer. If the remaining time until ERI falls under six months, an appropriate message is displayed.

Behavior in case of ERI

- In the dual-chamber modes, the pacemaker switches to single-chamber pacing when the elective replacement indication has been reached. This replacement mode varies according to the programmed pacing mode and is shown on the programmer.
- The following functions are deactivated when ERI has been reached:
 - Night program
 - Rate adaptation
 - Rate fading
 - Atrial overdrive pacing
 - 2:1 lock-in protection
 - PVC lock-in protection
 - PMT protection
 - Rate hysteresis
 - Lead check
 - Active capture control (ACC)
 - AV hysteresis
 - IEGM recordings
 - Statistics
- Pacing and sensitivity:

Magnet effect	Automatic
Magnet rate	80 ppm for 10 cycles directly after magnet application
Basic rate	Synchronous with basic rate reduced by 4.5 – 11% ^{a)}
Pulse width	Programmed values
Pulse amplitude	<ul style="list-style-type: none"> Programmed values with ACC deactivated (Off) If ACC was active (ON) before: last pacing threshold measured +1.2 V
Sensitivity	Programmed values

Magnet effect	Asynchronous
Magnet rate	Asynchronous with rate at 80 ppm
Basic rate	Asynchronous with rate at 80 ppm
Pulse width	Programmed values

Magnet effect	Asynchronous
Pulse amplitude	<ul style="list-style-type: none"> Programmed values with ACC deactivated (Off) If ACC was active (ON) before: last pacing threshold measured +1.2 V
Sensitivity	Programmed values

Magnet effect	Synchronous
Magnet rate	Synchronous with basic rate reduced by 4.5 – 11%
Basic rate	Synchronous with basic rate reduced by 4.5 – 11% ^{a)}
Pulse width	Programmed values
Pulse amplitude	<ul style="list-style-type: none"> Programmed values with ACC deactivated (Off) If ACC was active (ON) before: last pacing threshold measured +1.2 V
Sensitivity	Programmed values

^{a)} In the DDD(R), DDT(R), DDT(R)/A, DDT(R)/V, DDO(R), VDD(R), VDI(R), VDT(R), VVI(R), VVT(R), AAI(R), AAT(R), AOO(R) pacing modes, the pacing rate decreases by 11%. In the DDI(R), DDI/T(R), DVI(R), DVT(R) pacing modes, the VA interval increases by 11%, this reduces the pacing rate by 4.5 – 11%, depending on the programmed AV delay.

Expected service time after ERI

- The following tables show the mean and minimum values for the remaining service times between ERI (Elective Replacement Indication) and EOS (End of Service) for the standard program and for a program with high pulse energy.
- The information is based on a lead impedance of 500 Ohm, 100% pacing and the data of the battery manufacturer.
- For a lead impedance of 300 Ohm instead of 500 Ohm, these times decrease by max. 30%.

Philos II D, Philos II DR(-T) Expected service time (in months) in DDD(R) mode

Interval		Standard program	Program with high pulse energy
ERI to EOS	Mean value	9	8
ERI to EOS	Minimum value	7	6

Philos II SLR Expected service time (in months) in VDDR mode

Interval		Standard program	Program with high pulse energy
ERI to EOS	Mean value	8	8
ERI to EOS	Minimum value	6	6

Philos II S(R) Expected service time (in months) in AAI(R)/VVI(R) mode

Interval		Standard program	Program with high pulse energy
ERI to EOS	Mean value	8	8
ERI to EOS	Minimum value	6	6

Note: The expected service times can vary from those specified in the tables if the settings are different from those listed in the tables.

Explantation and Device Replacement

Explanted devices (pacemakers, defibrillators) can be sent to the local BIOTRONIK representative for proper, environmental friendly disposal.

⚠ CAUTION

Risk of infection due to biologically contaminated material

Explanted devices are biologically contaminated. Avoid infections using usual safety procedures.

⚠ CAUTION

Non-foreseeable risks of reuse

The device is only intended for one-time use. It must not be resterilized and reused.

Device replacement

If you replace the device and continue using leads that have already been implanted, please perform the following precautionary measure:

⚠ CAUTION

Technical failures due to damaged leads

While replacing a device check the leads that were already implanted with the previous device for integrity to ensure flawless functioning of the device system.

Returning the device

Proceed as follows:

Step	Action
1	Clean the device using a sodium hypochlorite solution with at least 1% chlorine.
2	Then thoroughly rinse the device.
3	Send the device to BIOTRONIK.

Cremation

Explant the device prior to cremation of a deceased pacemaker or ICD patient.

Cautionary Notes

Basic Information

The active device, the lead(s), and, if used, the lead extensions and adapters, become part of an artificial pacing system upon implantation. These components are responsible for the proper functioning of the artificial pacing system, which is based on the physiological condition of the patient.

⚠ CAUTION

The following notes are intended to emphasize some aspects that have been deemed especially important in the medical literature for evaluating and avoiding risks. This information could be useful in evaluating and avoiding risks, but it is not a substitute for the study of medical literature.

Notes for patients

Supplied with the device are a patient brochure and a patient ID card. The brochure is supplementary information for the patient and supports the instructions provided the physician. The patient brochure includes information on contraindicated methods of diagnosis and therapeutic procedures and also on special environmental conditions such as magnetic fields and electromagnetic interference fields.

Note: Give the patient's brochure and the completed patient ID card to the patient prior to discharge.

Possible Medical Complications

Overview

Possible medical complications of pacemaker therapy include the following:

- Formation of necrotic tissue
- Thrombosis, embolism
- Foreign body rejection phenomena
- Muscle and/or nerve stimulation
- Cardiac tamponade
- Infection
- Increase in pacing threshold
- Pacemaker-mediated arrhythmias

Actions to prevent device-induced complications

CAUTION

Possible unphysiological rhythm changes and asynchronous pacing upon magnet application

Please note that unphysiological rhythm changes can occur upon magnet application with programmed magnet mode asynchronous or auto. Asynchronous pacing occurs during magnet application for the programmed modes as follows:

- auto - 10 cycles for the duration of magnet application
- asynchronous - permanent for the duration of magnet application

CAUTION

Transmission of atrial tachycardia to the ventricle

Set the following parameters to avoid unphysiological pacing in the ventricle in cases of high atrial rates or sinus tachycardia:

- Activate Mode Switching for indicated patients.
- Set the upper rate and the refractory periods to prevent abrupt ventricular rate switching.
- Prefer Wenckebach response and avoid 2:1 behavior.

Set all parameters to prevent constant switching between atrial and ventricle controlled modes to prevent unphysiological rhythm switching upon loss of AV sequential pacing.

CAUTION

Pacemaker-mediated tachycardia in cases of retrograde conduction

Pacemaker-mediated tachycardia can occur in patients with retrograde conduction.

- Measure the retrograde conduction time.
- Switch on PMT protection to prevent pacemaker-mediated tachycardia.
- Set the VA criterion.

CAUTION

Unipolar pacing interferes when ICD is implanted

If an ICD is implanted at the same time as a pacemaker and a lead failure occurs, it is possible to switch to unipolar pacing after a pacemaker reset or using the automatic lead check function. This can falsely inhibit or trigger delivery of tachyarrhythmia therapy from the ICD.

- Switch off the lead check function.

Risky Therapeutic and Diagnostic Procedures

External defibrillation

The device is protected against the energy that is normally induced by external defibrillation. Nevertheless, any implanted device may be damaged by external defibrillation. Specifically, the current induced in the implanted leads may result in necrotic tissue formation close to the electrode/tissue interface. As a result, sensing properties and pacing thresholds may change.

Note: Place adhesive electrodes anterior-posterior or perpendicular to the axis formed by the device to the heart at least 10 cm away from the device and from implanted leads.

Due to possible damage to the patient or the device, the use of the following procedures is contraindicated:

- Therapeutic ultrasound and diathermy
 - Damage to the patient via excess warming of body tissue near the device system

- Radiation therapy
 - Shield the device sufficiently against radiation.
 - Check the system for integrity after applying radiation.
 - Radiation can cause latent damage.
- Transcutaneous Electrical Nerve Stimulation (TENS)
- Lithotripsy
- Magnetic resonance imaging and the associated magnetic flux density
 - Damage or destruction of the device system due to strong magnetic interaction
 - Damage to the patient via excess warming of body tissue near the device system
- Electrocautery and high-frequency surgery
 - Damage to the patient via the induction of arrhythmia or ventricular fibrillation
- Hyperbaric oxygen therapy
- Applied pressures higher than normal pressure

CAUTION

Interference with the device and risk to the patient due to electrical current during medical treatment

If electrical current from an external source is conducted through the body for diagnostic or therapeutic purposes, then the device has to be switched off or carefully monitored during the initial treatment phases.

Possible Technical Complications

Events that could compromise functioning include, for example:

- Pacemaker component error
- Battery depletion
- Lead dislocation
- Lead fracture
- Insulation defect

CAUTION

Technical failures due to damaged leads.

Switch the automatic lead check function on to detect lead errors. Impedance values that indicate a loss of lead integrity are documented in the event list.

Possible Interference

Possible impairments due to EMI

All devices (PM/ICD) are susceptible to interference. The signals generated by sources of interference could be interpreted as heart activity by the device, and/or could disturb measurements required for adapting the rate of the device:

- Depending on the pacing mode and the type of interference, these sources of interference may lead to pacemaker pulse inhibition or triggering, an increase in the sensor-dependent pacing rate, or a fixed-rate pulse delivery.
- Under unfavorable conditions, for example during diagnostic or therapeutic procedures, the interference sources may induce such a high level of energy into the artificial pacing system that the device and/or cardiac tissue around the lead tip is damaged.
- BIOTRONIK devices have been designed so that their susceptibility to EMI is minimized.
- However, due to the variety and intensity of EMI, absolute safety is not possible. It is generally assumed that EMI produces only minor symptoms, if any, in patients.

Note: If a device sensitivity value of < 2.0 mV/unipolar is set, interference from electromagnetic fields could result. Therefore, it is recommended that a value of ≥ 2.0 mV/unipolar be set according to paragraph 28.22.1 of the EN 45502-2-1 standard. Setting sensitivity values < 2.0 mV/unipolar requires explicit clinical need. Values like this can only be set and retained with physician supervision.

Note: Indicate the possibility of interference to your patient if this is expected to have a clinically relevant impact. Protect the device against the interference itself or its effect by means of appropriate programming.

Possible sources of interference

Interference can be caused by:

- Household appliances
- Safety locks or anti-theft installations
- Strong electromagnetic fields

- Cellular phones and transmitters (CardioMessenger)
 - Patients are advised to hold cellular phones to the ear opposite the side on which the device is implanted. Furthermore the cellular phone or transmitter should be held at least 15 cm away.
 - Cellular phones and transmitters emit signals when in stand-by mode, i.e. when not in use. Therefore these devices should not be carried in a breast pocket or within a radius of 15 cm around the device.
 - Electromagnetic interference only has a temporary effect. Devices function again properly when the cellular phone or transmitter is moved away.

Technical Data

Pacing Modes

Mode	D, DR(-T)	SLR	S, SR
DD(R)	x		
DDT(R)/A	x		
DDT(R)/V	x		
DDT(R)	x		
DDT(R)	x		
DDI(R)	x		
DDI(R)/T	x		
DVI(R)	x		
DVT(R)	x		
D00(R)	x		
VDD(R)	x	x	
VDT(R)	x	x	
VDI(R)	x	x	x
VVI(R)	x	x	x
VVT(R)	x	x	x
V00(R)	x	x	x
AAI(R)	x		x

Mode	D, DR(-T)	SLR	S, SR
AAT(R)	x		x
A00(R)	x		x
Off (temporary)	x	x	x

Pulse and Timing Parameters

- Basic rate: The corresponding intervals (t) correlate with the rates (r) according to the formula $t = 60,000 / r$ (t in ms, r in ppm).
- Rate hysteresis: not active during the night program.
- Rate adaptation: With Philos II DR(-T) the DDIR mode is only programmable in conjunction with activated mode switching.

Parameter	Value range	D, DR(-T)	SLR	S, SR
Basic rate	30 ... (1) ... 88 ... (2) ... 122 ... (3) ... 140 ... (5) ... 30 ... (1) ... 88 ... (2) ... 122 ... (3) ... 140 ... (5) ... 180 ppm	x	x	x
Night rate	Off; 30 ... (1) ... 88 ... (2) ... 122 ... (3) ... 140 ... (5) ... 180 ppm	x	x	x
Begin of night/end of night	00:00 ... (00:10) ... 00:00 ... (00:10) ... 23:50 h	x	x	x
Rate hysteresis	Off; -5 ... (-5) ... (-5) ... -80 ppm	x	x	x
Repetitive rate hysteresis	Off; 1 ... (1) ... 10 cycles	x	x	x
Scan rate hysteresis	Off; 1 ... (1) ... 10 cycles	x	x	x
Upper tracking rate (UTR)	100 ... (10) ... 100 ... (10) ... 140; 160; 185 ppm	x	x	x
Tachycardia mode	2:1; WKB (automatic setting)	x	x	
Rate limit	195 ... 220 ppm	x	x	x
IRsplus	Off; On	x	x	
Dynamic AV delay	Off, low; medium; high; individual; fixed	x	x	

Parameter	Value range	D, DRI (-T)	S, R	S, SR
AV delay	15; 50; 75; 100; 120 ... (10) ... 200; 225; 250; 300 ms (programmable in 5 rate ranges)	x	x	
AV hysteresis	Off; low; medium; high; IRSplus; negative	x	x	
AV repetitive hysteresis	Off; 1 ... (1) ... 6	x	x	
AV scan hysteresis	Off; 1 ... (1) ... 6	x	x	
Repetitive negative AV hysteresis	1 ... (1) ... 10 ... (3) ... 100 ... (10) ... 180	x	x	
Sense compensation	Off, -15 ... [-15] ... -120 ms	x	x	
AV safety delay	100 ms	x	x	
Cross channel blank A	56 ms	x	x	
Far field blanking (after Vs, Vp)	56 ... (25) ... 200 ms	x	x	
Cross channel blank V	16; 24; 32; 40; 48; 56; 72 ms	x		
Magnet effect	Auto; asynchronous; synchronous	x	x	x
Pulse amplitude A	0,1 ... (0,1) ... 4,8 ... (0,6) ... 8.4 V	x		
Pulse amplitude V	0,1 ... (0,1) ... 4,8 ... (0,2) ... 8.4 V	x	x	x ^{a)}
Pulse width A	0.1; 0.2; 0.3; 0.4; 0.5; 0.75; 1.0; 1.5 ms	x		
Pulse width V	0.1; 0.2; 0.3; 0.4; 0.5; 0.75; 1.0; 1.5 ms	x	x	x ^{a)}
Sensitivity A	0,1 ... (0,1) ... 1,5 ... (0,5) ... 7.5 mV	x	x	
	0,4 ... (0,4) ... 6.0 mV			x ^{a)}
Sensitivity V	0,5 ... (0,5) ... 7.5 mV	x	x	x
Refractory period A	200 ... (25) ... 775 ms	x	x	
Refractory period V	170; 195; 220; 250; 300; 350; 400 ms	x	x	x ^{a)}
Atrial refractory period extension	0 ... (50) ... 350 ms	x	x	x
Max. activity rate	80 ... (5) ... 180 ppm	x	x	x
Sensor gain	Auto; 1 ... 40 (in 32 steps)	x	x	x

Parameter	Value range	D, DRI (-T)	S, R	S, SR
Automatic sensor gain	Off; On	x	x	x
Sensor threshold	Very low; low; medium; high; very high	x	x	x
Rate increase	1; 2; 4; 8 ppm/cycle	x	x	x
Rate decrease	0.1; 0.2; 0.5; 1.0 ppm/cycle	x	x	x
Tachycardia behavior	Off, mode conversion, mode switching	x	x	
Mode conversion	Off; on [in modes DDD(R), DDT(R)/A, DDT(R)/V, and VDD(R)]	x		
	Off; on [in mode VDD(R)]		x	
Mode switching	Off; on [in modes DDD(R), DDT(R)/A, DDT(R)/V, and VDD(R)]	x		
	Off; on [in mode VDD(R)]		x	
Intervention rate	110 ... (10) ... 250 ppm	x	x	
Mode switching: Activation criterion	3 ... (1) ... 8	x	x	
Mode switching: Deactivation criterion	3 ... (1) ... 8	x	x	
PMT protection	Off; On	x	x	
VA criterion	250 ... (10) ... 500 ms	x	x	
Min. PVARP	Off; On	x	x	
Polarity sense A	Unipolar; bipolar	x		
	Bipolar		x	
Polarity sense V	Unipolar; bipolar	x	x	x
Polarity pace A	Unipolar; bipolar	x		
Polarity pace V	Unipolar; bipolar	x	x	x
Lead check	Off; On	x	x	x
ACC	Off; On	x	x	x
ACC: min. amplitude	0.1 ... (0.1) ... 4.8 ... (0.2) ... 6.4 V	x	x	x
ACC: max. amplitude	2.4; 3.6; 4.8; 6.4 V	x	x	x
Threshold search	Interval: 0.1; 0.3; 1; 3; 6; 12; 24 h; times of day: 0:00 to 24:00 h	x	x	x

Parameter	Value range	D, DR(-T)	SLR	S, SR
Safety	0.3 ... [0.1] ... 0.5 ... [0.1] ... 1.2 V	x	x	x
IEGM recordings	12 recordings, each max. 10 s	x	x	x
IEGM trigger	Mode switching or high atrial rate	x	x	
	High ventricular rate			x

a) Value range applies to A, if the atrial mode is set.

Preset Programs

Parameter/function	Programs:			D, DR(-T)	SLR	S, SR
	Factory	Standard	Safe			
Mode (DR(-T), D)	DDD	DDD(R) ^{a)}	VI	x		
Mode (SLR)	VDD	VDD	VI		x	
Mode (SR, S)	VVI	VVI(R) ^{b)}	VVI ^{c)}			x
Basic rate	60 ppm	60 ppm	70 ppm	x	x	x
Night program	Off	Off	Off	x	x	x
Rate hysteresis	Off	Off	Off	x	x	x
Upper rate	130 ppm	130 ppm	—	x	x	
IRSpplus	Off	Off	—	x	x	
Dynamic AV delay	Low	Low	—	x	x	
AV hysteresis	Off	Off	—	x	x	
Sense compensation	-45 ms	-45 ms	—	x	x	
AV safety delay	100 ms	100 ms	—	x	x	
Cross channel blank A	56 ms	56 ms	—	x	x	
Cross channel blank V	32 ms	32 ms	—	x		
Magnet effect	Auto	Auto	Auto	x	x	x
Pulse amplitude A	3.6 V	3.6 V	—	x	x	
Pulse amplitude V	3.6 V	3.6 V	4.8 V	x	x	x

Parameter/function	Programs:			D, DR(-T)	SLR	S, SR
	Factory	Standard	Safe			
Pulse width A	0.4ms	0.4ms	—	x		
Pulse width V	0.4ms	0.4ms	1.0ms	x	x	x
Sensitivity A	1.0 mV	1.0 mV	—	x	x	
Sensitivity V	2.5 mV	2.5 mV	2.5 mV	x	x	x
Refractory period A	425 ms	425 ms	—	x	x	
Refractory period V	250 ms	250 ms	300 ms	x	x	x
ARP extension	0 ms	0 ms	—	x	x	
Mode conversion	Off	Off	—	x	x	
Mode switching	Off	On	—	x	x	
Mode switching: activation criterion	—	5-out-of-8	—	x	x	
Mode switching: deactivation criterion	—	5-out-of-8	—	x	x	
Intervention rate	—	160	—	x	x	
Far field blanking	56 ms	56 ms	—	x	x	
Switch to	—	DDIR	—	x	x	
Min. PVARP	—	235 ms	—	x	x	
Sensor threshold	—	Medium	—	x	x	
Sensor gain	—	4	—	x	x	
Automatic sensor gain	—	On	—	x	x	
Rate increase	—	2 ppm/s	—	x	x	
Max. activity rate	—	120 ppm	—	x	x	
Rate decrease	—	0.5 ppm/s	—	x	x	
Polarity pace A	Unipolar	Unipolar	Unipolar	x		
Polarity pace V	Unipolar	Unipolar	Unipolar	x	x	x
Polarity sense A	Unipolar	Unipolar	Unipolar	x		
Polarity sense V	Unipolar	Unipolar	Unipolar	x	x	x
PMT protection	Off	On	—	x	x	

Parameter/function	Programs:			D, DR(-T)	SLR	S, SR
	Factory	Standard	Safe			
VA criterion	—	380 ms	—	x	x	
Lead check	Off	On	Off	x	x	x
Auto-initialization	On	—	—	x	x	x
ACC	Off	Off	—	x	x	
IEGM recording	Off	Off	Off	x	x	

a) After auto-initialization, the device is set to DDD mode.

b) After auto-initialization, the device is set to VVI mode.

c) When the device is set to AAI mode, the AAI safe program is used.

Material in Contact with Human Tissue

Component	Material
Housing	Titanium
Sealing plug	Silicone
Header	Epoxy resin
Coating (if used)	Silicone

Programmable

Programmable that can be used:


Programmable
PMS 1000plus
PMS 1000 C
PRT 1000
TMS 1000plus
TMS 1000
ICS 3000 Implant Control System

Electrical Characteristics

Electrical characteristics*)

Circuit	Hybrid electronics with VLSI-CMOS chip
Input impedance A	>10 kOhm
Input impedance V	>10 kOhm
Pulse form	Biphasic, asymmetric
Polarity	Cathodic

Pulse waveform graph:

The pulse amplitude reaches its maximum value at the beginning of the pulse (U_a). With increasing pulse width (t_b), the pulse amplitude is reduced dependent on the pacing impedance.

Power consumption	DR	D	SLR	SR	S
BOS, inhibited	13 μ A	13 μ A	13 μ A	13 μ A	13 μ A
BOS, 100% pacing	21 μ A	21 μ A	16.5 μ A	16.5 μ A	16.5 μ A

Electrically conductive housing surface	
Uncoated	32.8 cm ²
Coated	7.23 cm ²
Electrically conductive housing shape	
Uncoated	Flattened ellipsoid
Coated	Ellipse

*) 37°C, 500 Ohm

Battery

Type	Li/I	
Manufacturer	Wilson Greatbatch	Litronik
Model	WG 8431	LIS 3150
Open-circuit voltage	2.8 V	
Nominal capacity ^{a)}	1.3 Ah	

*] Battery manufacturer data

Service Times

Service time ^{a)} in years	DR(-T)	D	SLR	SR	S
Nominal service time ^{b)} at pulse amplitudes 3.6 V	7.2	7.2	10.0	10.0	10.0
Expected service time ^{c)} at pulse amplitudes 3.6 V	5.9	5.9	7.3	7.3	7.3
Remaining capacity at ERI (in Ah)	0.1	0.1	0.1	0.1	0.1

a) Device set to standard program, 100% pacing, calculated with the battery manufacturer data

b) Calculated according to the formula: $T = 2740 \times C_{bat} / (IBOS + IEOS)$

c) Expected service time based on all available data

Mechanical Characteristics

	DR, D, SLR	DR-T	SR, S
Weight [g]	26	27	25
Volume [cm ³]	12	12	11
Dimension [mm]	6 x 43 x 53	6 x 44 x 51	6 x 39 x 53

X-ray identification

Implant	
Philos II: DR, D, SLR, SR, S	ET
Philos II: DR-T	KP

Storage Conditions

Values	Value range
Relative humidity	Max. 70%
Temperature	5–55°C
Pressure	0.7–1.5 bar

Estimated Tolerances of the Factory Program

- Factory settings at 37°C, 500 Ohm
- Data according to EN 50061

D(R/-T), SLR, S(R)

Basic rate, interference rate	60 +/- 1.5 min ⁻¹
Basic interval	1000 +/- 20 ms
Escape interval	1000 +/- 20 ms
Magnet rate	90 +/- 3 min ⁻¹ (for 10 cycles)
Magnet interval	664 +/- 20 ms (for 10 cycles)
AV delay	
Basic rate	180 +/- 15/-5 ms
< 70 ppm	180 +/- 15/-5 ms
70–90 ppm	160 +/- 15/-5 ms
91–110 ppm	140 +/- 15/-5 ms
111–130 ppm	120 +/- 15/-5 ms
> 130 ppm	100 +/- 15/-5 ms

	Atrium	Ventricle
Pulse amplitude	3.6 + 0.1/-0.7	3.6 + 0.1/-0.7
Peak value		
Pulse amplitude	3.3 + 0.1/-0.7	3.3 + 0.1/-0.7
EN 50061: Mean value		
Pulse width	0.42 +/- 0.02 ms	0.42 +/- 0.02 ms
Refractory period ^{a)}	425 + 10/-20 ms	250 + 10/-20 ms
Rate limit ^{a)}	200 + 20/-5 min ⁻¹	200 + 20/-5 min ⁻¹

^{a)} Does not apply to S(R): see S(R), p. 18.

D(R/-T)

Sensitivity	Atrium	Ventricle
15 ms sin ²	1.0 +/- 0.5 mV	-
40 ms sin ²	-	2.5 +/- 0.5 mV
Delta pulse	1.8 +/- 0.5 mV	2.5 +/- 0.5 mV

SLR

Sensitivity	Atrium	Ventricle
15 ms sin ²	0.2 +0.05/-0.1 mV	-
40 ms sin ²	-	2.5 +/- 0.5 mV
Delta pulse	0.24 +0.05/-0.1 mV	2.5 +/- 0.5 mV

S(R)

Sensitivity	Atrium	Ventricle
15 ms sin ²	2.0 +/- 0.5 mV	-
40 ms sin ²	-	2.5 +/- 0.5 mV
Delta pulse	2.5 +/- 0.5 mV	

Refractory period	300 + 10/-20 ms
Rate limit	200 + 20/-5 min ⁻¹

Product Line

Model	Lead connection	Order number
Philos II DR		
Uncoated	IS-1	341 826
Coated	IS-1	341 821
Philos II DR-T		
Uncoated	IS-1	343 175
Coated	IS-1	343 176
Philos II D		
Uncoated	IS-1	341 825
Coated	IS-1	341 820
Philos II SLR		
Uncoated	IS-1	341 822
Coated	IS-1	341 816
Philos II SR		
Uncoated	IS-1	341 824
Coated	IS-1	341 815
Philos II S		
Uncoated	IS-1	341 823
Coated	IS-1	341 819

Legend for the Label

The label icons symbolize the following:

Symbol	Meaning
	Manufacturing date
	Use by date, only valid if sterile package is undamaged! Never use outdated devices!
	Storage temperature
REF	BIOTRONIK order number
SN	Serial number
PID	Product identification number
STERILE EO	Sterilized with ethylene oxide
	Do not resterilize!
	Single use only. Do not reuse!
	Non-sterile
	Consult operating instructions!
	Contents

Symbol	Meaning
	Do not use if packaging is damaged!
CE	European approval mark
	Transmitter with non-ionizing radiation on specified frequency; turned off at delivery
 5/6 mm	Upsizing adapter sleeve for leads with PIN-Lock PE (5 mm) connector for connection to pacemakers with PEC (6 mm) connector
 IS-1 UNI/BI	Bradycardia therapy device with NBG Code and designation of compatible leads (example)
 IS-1 UNI/BI	Silicone-coated bradycardia therapy device with NBG Code and designation of compatible leads (example)
	Screwdriver
UNI / BI	Unipolar/bipolar polarity configuration
Vp; Vs	Ventricular paced event; Ventricular sensed event
Ap; As	Atrial paced event; Atrial sensed event

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All rights reserved. Specifications subject
registered trademarks of BIOTRONIK SE & Co. KG