

Nominal Specifications		
	1888TC (Straight)	1882TC (Atrial J)
<b>Connector type:</b>	IS-1 <sup>1</sup> bipolar	
<b>Lead introducer (minimum recommended)</b>	6 F <sup>2</sup> (without retained guidewire)	
<b>Maximum diameter of lead body:</b>	1.87 mm	2.00 mm
<b>Fixation mechanism:</b>	helix	
<b>Electrode configuration:</b> Tip Ring	active helix and active mapping collar cylindrical	
<b>Electrode spacing:</b> Tip to ring	10 mm	
<b>Electrode surface area:</b> Tip Ring	8.5 mm <sup>2</sup> 16 mm <sup>2</sup>	
<b>Electrode length:</b> Tip (extended)	2 mm	
<b>Electrical resistance:</b> Tip to connector pin Ring to connector ring	19 Ω 39 Ω	19 Ω 24 Ω
<b>Materials:</b> Conductors Connectors Lead body Suture sleeve Tip and ring electrodes Mapping collar Steroid-eluting plug  Stylet	multifilar MP35N <sup>3</sup> stainless steel Optim <sup>®</sup> 4 insulation with exterior Fast-Pass <sup>®</sup> coating silicone rubber titanium nitride-coated platinum-iridium alloy platinum-iridium alloy silicone rubber with less than 1.0 mg dexamethasone sodium phosphate stainless steel	

Nominal Specifications		
	1888TC (Straight)	1882TC (Atrial J)
<b>Typical number of rotations to extend helix for initial placement<sup>5</sup>:</b>		
Straight stylet	6 - 11 <sup>6</sup>	9 - 14 <sup>6</sup>
J stylet	9 - 14 <sup>6</sup>	
<b>Stylet diameter (of those included in package):</b> X-soft (light green knob) Soft (green knob) J-shaped soft (green knob/ white button)	0.35 mm (0.014 in)	

1. St. Jude Medical IS-1 lead connectors are compatible with connector cavities that conform to the international connector standard ISO 5841-3.
2. Compatible with introducers such as those manufactured by Daig or Pressure Products.
3. MP35N is a registered trademark of Standard Pressed Steel Technologies.
4. A silicone-polyurethane copolymer.
5. Lead length, stylet configuration, and differences in anatomy may cause variation in the number of rotations required to extend the helix.
6. Not to exceed 35 rotations.

Accessories
Stylet Kit with Clip-On Tool
Stylet Kit with Fixation Tool
Locator <sup>®</sup> Steerable Stylet <sup>1</sup>

1. Model 1888TC only.

Cardiac Rhythm Management Division

[sjm.com](http://sjm.com)



For indications, contraindications, warnings, and cautions, see the accompanying manual.



November 2012  
Art 60008638/B